

 **ROGG**<sup>®</sup>vet

Product catalogue 2022/2023

 + 49 (0) 8133/9320-10  
 + 49 (0) 8133/9320-25  
 [sales@roggvet.de](mailto:sales@roggvet.de)  
 [www.roggvet.de](http://www.roggvet.de)





Innovative  
products.



Strong brands.



Sustainable  
animal health.

## On behalf of humans and animals

For more than 50 years, ROGG has been a pioneer in modern wound treatment in the field of human medicine and, for the past five years, also in veterinary medicine with our quality brand EQUIMOLL®, and our partner brand benartis.

With ROGGvet – our young veterinary division – we are guided by the same credo for animals as for people: Innovation and health are inextricably linked.

At ROGGvet, we have made it our mission to provide you with contemporary, economical solutions so that you could meet all the challenges of everyday life, the demands of social change as well as the needs of your animal patients.

Reliability, innovation, and sustainability are our guiding principles that define our services, our product selection, and our internal operations.

## Everything from a single source: wound treatment, hygiene, practice equipment

On the basis of many years of experience and in direct cooperation with you, our customers, and our extensive network of specialists, we are constantly improving our products and brands – always looking for more innovative brands and economic solutions with the aim of matching your needs even better and our own requirements in terms of progressiveness and sustainability.

In our current ROGGvet catalogue you will find a product range of strong brands matching the human medicine quality level and combining economic efficiency and innovation. Our portfolio, which is tailored both to large and/or cloven-hoofed and small animals, covers not only wound treatment and fixation with bandages and dressings, but also the areas of hygiene and practice demand.

 ROGG<sup>®</sup>vet



## Your reliable partner

At ROGGvet, we don't see ourselves just as a retailer of veterinary products, but as your reliable partner for animal health. Our customer service and our sales team know your everyday challenges and are pleased to be at your disposal for questions, personal advice, and suggestions.

We deliver the products you need within a few days, and in case something doesn't go according to plan, we promise you open communication and quick and fair solutions.


Fairness and transparency is not just a principle that guides us in our customer relationships. ROGGvet is a division of ROGG Verbandstoffe GmbH & Co. KG. Everyone working here is part of a diverse team that offers equal opportunities, family-friendliness, and individual work models – regardless of nationality, conviction, origin, skin color, gender, age, lifestyle, and everything else that makes people who they are.

## Project "Sustainable Animal Practice"

Environmental protection and sustainability are considered major challenges for medicine, especially under the steadily increasing economic pressure – also in veterinary practices. We, at ROGGvet, believe that responsible use of resources, economic efficiency, and modern medical care can not only be combined but must be brought together. In concrete terms, this means: shortening treatment times through particularly innovative approaches, reducing packaging waste, and minimizing the ecological footprint. Because humans and animals are dependent on an intact environment.

ROGGvet contributes to the "Sustainable Veterinary Practice" project. We are constantly optimizing our own workflows and processes with the aim of optimal future viability. We ship climate-neutrally with GLS and get 100 percent green electricity from "Lichtblick". Last year, we saved 8.72t of CO<sub>2</sub> through the purchase of green energy alone. And due to further optimization, we will soon be able to reduce our packaging waste by 80 percent. And we won't stop there.

Any product missing from our portfolio? Please give us a call or write an email! We would be happy to order the items you need for you.







## Bandages and dressings

### Gauze rolls 14

EQUIMOLL® ORIGINAL	14
EQUIMOLL® ORIGINAL Sample	15
EQUIMOLL® JUNIOR	15
benartis cotton wool roll Dainamic	16

### Compresses 16

ROGG gauze compresses	16
BeeSana gauze compresses	17
Zemuko® cellulose non-woven compress	17
ROGG gauze swabs	18

### Fixation bandages 18

ROGGELAST	18
EQUILAST	19
EQUIHAFT	19
ROGGHAFT Crepe Color	20
ROGGHAFT Color	21
benartis F Dainamic	22
Valueline Flex Bandage	24

### Zinc paste bandages 24

EQUIZINK	24
----------	----

### Padding material 26

EQUI-SOFT	26
benartis cotton wool Dainamic 50/50	26
ROGG synthetic cotton bandage	27
benartis cotton wool Xactic	27
ROGG foam bandage	28
snögg AnimalPOLSTER	28

### Adhesive bandages 30

EQUI-PLAST	30
benartis A Plaster Bandage	30

### Plasters 31

ROGG Stretch	31
Miro Plast	32
ROGG Silk	32

### Tubular bandages 34

EQUICOMP	34
ROGG terry cloth tubes	34



## Wound management

### Wound cleansing 36

KerraSol™ wound irrigation solution	36
KerraSol™ wound spray gel	36

### Wound dressings 37

ROGG Bio-Kollagen (organic collagen)	37
ROGG Vulcosan AMC	37
curea clean Veterinary	38
Revamil® wound dressing	38

### Ointments 39

Revamil® gel	39
Revamil® balm	39
NAWAlution Gel	40
ROGG wound & burn ointment	40

### Wound foam dressings 41

snögg AnimalVET	41
snögg AnimalSOFT	41



## Disinfection & hygiene

### Skin and hands 42

ROGG washing lotion	42
ROGG skin/hand disinfection	42
ROGG hand disinfectant gel	43

### Skin, hands and surfaces 44

Isopropyl alcohol 70 %	44
Isopropyl alcohol > 99 %	44

### Skin – preoperative 45

kodan® tincture forte, colorless	45
kodan® tincture forte, colored	45

### Surfaces 46

ROGG disinfectant wipes AF Q*	46
ROGG rapid disinfection set AF plus**	46
ROGG system non-woven wipes	47
ROGG rapid disinfection AF plus	48

### ROGG rapid disinfection AF Q 48

ROGG surface disinfection AF	49
ROGG surface disinfection AF Q	49

### Instruments 50

ROGG instrument disinfection AF	50
ROGG instrument disinfection VV	50

### Accessoires 51

Disinfection tray	51
Application aids	51
Dosing bottle	51
Hand spray head for spray disinfection	51
Dosing pump	51
Canister pump	51
Canister tap	51



## Practice & OR supplies

### Protective clothing 54

AMPri Style by Med Comfort	55
PE veterinary gloves AMPri Med-Comfort	56
AMPri Eco-Plus Latex	56
AMPri Blue Eco-Plus Vinyl	57
Medi-Inn Nitril	57
Gentle Skin® Premium OP™	58
Reference OP®	58
Suavel® Protec Plus	59
Suavel® Comfort Plus	59
Suavel® Antifog	60

### Reusable gown "Krefeld" 61

BeeSana® PP protective sheath	61
BeeSana® PP gown	62
BeeSana® SMS gown	62
Suavel® Astrid	63
Surgical overshoes	63

### Skin care 64

VeraPanthen™	64
ROGG Skin Protection Cream	64

### Accessoires 65

### Fripa towel dispenser 65

Fripa disposable towels	65
Ecotex Plus disposable towels	65
ROGG dispenser box	65

### Adhesive tapes 66

Certoplast Prima insulating tape	66
tesa textile adhesive tape	66
tesa masking tape	67
tesa extra Power® Universal	67
benartis silver tape professional	68
benartis duct tape	68

### Surgical drapes 69

### mediware sterile surgical drapes 69

ZARYS incision drape	69
----------------------	----

### Underlays 70

Pulp / cellulose	70
Lohmann pulp	70
curea SafeCare	71



## Care

### Skin and Fur 72

DERMAVET LAVAPIROX® shampoo	72
DERMAVET LAVAPIROX® foam	73
DERMAVET LAVAPIROX® fluid	73
DERMAVET LAVAPIROX® cool gel	74

### Ears and Teeth 75

DERMAVET LAVAPIROX® clean	75
DERMAVET LAVAPIROX® vet anti-plaque C	75

### Feed Supplements 76

Equital Arthro Star Plus Liquid	76
Equital Broncho Horse Liquid	76
Equital Bone Star Powder	77

### Equital Hepa Forte Extra 78

Equital Dermasan Horse Liquid	78
Equital Immun Aktiv Liquid	79





Our  
own brand



EQUIMOLL® – a brand that stands for tradition, competence and premium quality in the field of horse injury treatment. EQUIMOLL® has been part of our ROGG Health Innovation brand family since 2017. Since then, we have been able to steadily expand our range around our best-known product EQUIMOLL® ORIGINAL by developing high-quality products in the area of bandages and dressings.

The highest quality standards demanded by our partners and our EQUIMOLL® products guarantee tried and tested developments, responsible, sustainable use of materials and consistently first-class quality.

# benartis

In 2019, the young brand benartis joined the portfolio of our partner company, Hemel Vertriebs-GmbH, in the veterinary sector and is marketed by ROGG. benartis combines our demands on quality and effectiveness with groundbreaking, trend-oriented functionality and an excellent price-performance ratio. The health and well-being of the animals are our top priority. At the same time, smart, intuitive materials ensure convenient and easy handling.



Our  
partner brand





# The perfect solution for every problem

## Application examples



### SUPERFICIAL WOUNDS



**Zemuko® cellulose non-woven compress**

**MORE ON PAGE 17**



**ROGG gauze compress**

**MORE ON PAGE 16**



**ROGG Stretch**

**MORE ON PAGE 26**

### MODE OF ACTION

- Protection against contamination
- Promotion of the healing process

### ADDITIONAL PRODUCT RECOMMENDATIONS

for better wound healing:

- NAWAlution Gel
- Revamil® gel
- Revamil® balm

### CANNON BONE & JOINT BANDAGE



**EQUIMOLL® ORIGINAL**

**MORE ON PAGE 14**



**benartis F Dainamic**

**MORE ON PAGE 22**

### MODE OF ACTION


- as a support bandage, for stability
- as a sprue bandage
- as a wet-to-dry dressing with a moist wound pad

### ADDITIONAL PRODUCT RECOMMENDATIONS

for a wet-to-dry dressing:

- ROGG Vulcosan AMC
- ROGG Bio-Kollagen
- curea clean
- Revamil® wound compress

### SCURF, THRUSH



**tesa textile adhesive tape**

**MORE ON PAGE 66**



**EQUIMOLL® ORIGINAL**

**MORE ON PAGE 14**



**benartis F Dainamic**

**MORE ON PAGE 22**

### MODE OF ACTION

- Protection against disease progression
- Protection against ingress of moisture, bacteria and dirt
- Promotion of the healing process



## The perfect solution for every problem

### Application examples



#### CANNON BONE & JOINT BANDAGE



**EQUIMOLL® ORIGINAL**

**MORE ON PAGE 14**



**EQUILAST**

**MORE ON PAGE 19**



**tesa masking tape**

**MORE ON PAGE 67**

#### MODE OF ACTION

- applicable on the horse's leg in a variety of ways

#### JOINT BANDAGE



**EQUILAST**

**MORE ON PAGE 19**



**tesa masking tape**

**MORE ON PAGE 67**



**EQUI-PLAST**

**MORE ON PAGE 30**

#### MODE OF ACTION

- as a support bandage, for stability
- as a sprue bandage
- as a wet-to-dry dressing with a moist wound pad

#### JOINT BANDAGE



**EQUIHAFT**

**MORE ON PAGE 19**

#### MODE OF ACTION

- applicable on the horse's leg and other parts of the body in a variety of ways for fixation





Our most popular item



## EQUIMOLL® ORIGINAL

Handmade cotton wool compress roll

EQUIMOLL® ORIGINAL is a pure cotton wool bandage, handmade in Germany, made from a premium cotton-viscose mix, wrapped in gauze on one side. Fluffy – soft with an uniform thickness, does not stick to the wound and does not leave any fluff when cut to size. Optimal width for use with large animals.

Ingredients: cotton, viscose

size	colour	pcs./pack. unit	Item No.
40 cm x 6 m	white	1 roll	203589

Carton units: 8 rolls, 192 rolls (1 pallet)

### ADVANTAGES

- lint-free wound covering thanks to the interwoven gauze structure
- high economic efficiency, since it can be used multiple times and can be cut exactly to size
- no "fluffing"
- shortened, sustainable treatment duration
- low material consumption
- high softness thanks to high quality viscose
- gauze is utilized on one side and not in the middle
- high stability
- no wrinkling / no lacing
- extra long – 6 m

### APPLICATION AREAS

- ideal for hoof and sprue bandages
- ideal for protecting injuries, bandaging and padding open and large wounds
- for creating cast bandages
- to keep tendons or muscles warm during training or transport
- head bandages
- tubular bandages of any kind
- for padding support bandages (e.g. Robert Jones bandages)
- tail regions dressings



Available in test size

And for the little ones



## EQUIMOLL® ORIGINAL Sample

The original in test size

size	colour	pcs./pack. unit	Item No.
40 cm x 3 m	white	1 roll	204699

Carton units: 8 rolls

ONLY WHILE SUPPLIES LAST!

## EQUIMOLL® JUNIOR

For ponies, foals and cloven-hoofed animals

size	colour	pcs./pack. unit	Item No.
20 cm x 6 m	white	2 rolls	203718

Carton units: 16 rolls

ONLY WHILE SUPPLIES LAST!



**Supplementary product:**  
EQUIHAFT, benartis F Dainamic





benartis cotton wool roll Dainamic

Compress roll

This gauze roll or compress roll is a cotton wool bandage made of pure cotton, wrapped in gauze. It can be cut to size and does not stick to the wound. 950g at the specified size.

Ingredients: 100 % cotton

size	colour	pcs./pack. unit	Item No.
40 cm x 5 m	white	1 piece	1109

Carton units: 10 rolls, 120 rolls

APPLICATION AREAS

- for protection against injuries
- for bandaging and padding open and large wounds
- for keeping tendons or muscles warm during training or transport



Supplementary product:  
benartis F Dainamic



ROGG gauze compresses

Non-sterile, 8x folded, 17 threads

Gauze compresses with cut edges folded in and good absorbency. Soft, air-permeable wound pad made of pure cotton, 17 threads, 8 layers.

Ingredients: 100 % cotton

size	colour	pcs./pack. unit	Item No.
5 cm x 5 cm	white	100 pieces	731
7.5 cm x 7.5 cm	white	100 pieces	735
10 cm x 10 cm	white	100 pieces	734

Carton units: 1,000 pieces

APPLICATION AREAS

- for the care of wounds of all kinds
- as a wound pad



Supplementary product:  
EQUILAST, ROGGEAST



BeeSana gauze compresses

Non-sterile, 12x folded, 17 threads

Gauze compresses with cut edges folded in and strong absorbency. Soft, air-permeable wound pad made of pure cotton, 17 threads, 12 layers.

Ingredients: 100 % cotton

size	colour	pcs./pack. unit	Item No.
5 cm x 5 cm	white	100 pieces	732
7.5 cm x 7.5 cm	white	100 pieces	736
10 cm x 10 cm	white	100 pieces	727

Carton units: 1,000 pieces

size	colour	pcs./pack. unit	Item No.
10 cm x 20 cm	white	100 pieces	780

Carton units: 500 pieces

APPLICATION AREAS

- for the care of wounds of all kinds
- as a wound cover
- as a wound pad

ADVANTAGES

- particularly absorbent



Supplementary product:  
EQUILAST, ROGGEAST



Zemuko® cellulose non-woven compress

Very soft, with one-sided, non-sticking fleece coating

Soft, economical universal compress, consisting of a fine fleece coating with an absorbent layer of highly bleached cellulose. The wound-friendly, non-woven fabric and the absorbent layer are connected to one another without creases by dot grids. It can therefore be cut from the roll firmly and lint-free to edge.

Ingredients: cellulose, polyester, viscose, dispersion adhesive

size	colour	pcs./pack. unit	Item No.
10 cm x 10 m	white	1 piece	774

Carton units: 20 rolls

APPLICATION AREAS

- for the absorption of exudate
- as an additional absorbent layer
- for ointment dressings and moist compresses
- for padding
- for laying out instrument trays



Supplementary product:  
ROGG Dispenser Box  
ROGG Stretch





ROGG gauze swabs

Highly absorbent swabs

Swabs made of gauze, pure cotton, white, 20 threads

Ingredients: 100 % cotton

size	colour	pcs./pack. unit	Item No.
Size 1 - plum-sized	white	2x 250 pieces	1131
Size 2 - hazelnut-sized	white	500 pieces	1132
Size 3 - walnut-sized	white	500 pieces	11321
Size 4 - egg-sized	white	500 pieces	11322

Carton units: 1,000 pieces

APPLICATION AREAS

- for wound care



ROGELAST

Soft, elastic gauze bandage

This edged, white and non-slip fixation bandage can be applied easily, wrinkle-free and in a material-saving manner. It is approx. 4 m long (stretched).

Ingredients: 40 % viscose, 60 % polyamide

size	colour	pcs./pack. unit	Item No.
10 cm x 4 m	white	20 pieces	1449

Carton units: 240 pieces

size	colour	pcs./pack. unit	Item No.
12 cm x 4 m	white	20 pieces	1450

Carton units: 440 pieces

APPLICATION AREAS

- for wrinkle-free fixation of wound dressings

ADVANTAGES

- non-slip and wrinkle-free



Supplementary product:  
ROGG Silk, tesa masking tape



EQUILAST

Elastic ideal bandage

EQUILAST is a white, bleached, universal bandage made from a cotton-polyamide-elastane blend. It has an elasticity of approx. 90 % and can absorb small amounts of sweat.

Ingredients: 64 % cotton, 34 % polyamide, 2 % elastane

size	colour	pcs./pack. unit	Item No.
10 cm x 5 m stretched	white	12 pieces / box	204262

Carton units: 87 pieces

size	colour	pcs./pack. unit	Item No.
12 cm x 5 m stretched	white	8 pieces / box	204263

Carton units: 72 pieces

size	colour	pcs./pack. unit	Item No.
15 cm x 5 m stretched	white	10 pieces / box	204264

Carton units: 54 pieces

APPLICATION AREAS


- for a wrinkle-free application of bandages
- for fixing wound dressings on the horse's neck, head or leg

ADVANTAGES

- it is very comfortable to wear thanks to the fine, soft structure and provides optimal fixation



Supplementary product:  
ROGG Silk, tesa masking tape



EQUIHAFT

Elastic, self-adhesive fixation bandage

EQUIHAFT is an off-white fixation bandage, made of viscose thread and polyamide, coated with natural latex at certain points. It has an elasticity of approx. 75 %, is non-slip, air-permeable and does not stick to the skin or hair of the animals.

Ingredients: viscose, polyamide, latex

size	colour	pcs./pack. unit	Item No.
6 cm x 4 m	white	100 Pieces	203592

6 cm x 4 m white 10 Pieces 204265

8 cm x 4 m white 100 Pieces 203593

8 cm x 4 m white 10 Pieces 204266

10 cm x 4 m white 100 Pieces 203594

10 cm x 4 m white 10 Pieces 204267

12 cm x 4 m white 100 Pieces 203595

12 cm x 4 m white 8 Pieces 204268

12 cm x 20 m white 34 Pieces 203596

12 cm x 20 m white 2 Pieces 204270

APPLICATION AREAS

- for fixing wound dressings and padding material
- for applying compression and support bandages especially to moving parts of the body

ADVANTAGES

- particularly good adhesion and yet air-permeability thanks to optimal natural latex processing / distribution
- natural latex does not stick to skin or hair and, thus, enables an atraumatic dressing change
- enables high compression thanks to its 75 percent elasticity



Supplementary product:  
EQUIMOLL® ORIGINAL





ROGGHAFT Crepe Color

Cohesive, elastic fixation bandage with creped structure

This self-adhesive, non-slip, air-permeable and extremely skin-friendly fixation bandage can be gas and steam sterilized and has a creped structure that ensures optimum stability. It can be stretched to about 85 %. Stretched, it is about 20 m long. Ideal for the fixation of wound dressings and for the post-treatment of distortions.

Ingredients: 39 % cotton, 61 % polyamide, with latex coating

size	colour	pcs./pack. unit	Item No.
6 cm x 20 m	red	1 piece	14811
6 cm x 20 m	dark blue	1 piece	14812
6 cm x 20 m	green	1 piece	14813
6 cm x 20 m	purple	1 piece	14814
6 cm x 20 m	yellow	1 piece	14815
6 cm x 20 m	royal blue	1 piece	14816
6 cm x 20 m	royal black	1 piece	14817

Carton units: 16 pieces

size	colour	pcs./pack. unit	Item No.
8 cm x 20 m	red	1 piece	14821
8 cm x 20 m	dark blue	1 piece	14822
8 cm x 20 m	green	1 piece	14823
8 cm x 20 m	purple	1 piece	14824
8 cm x 20 m	yellow	1 piece	14825
8 cm x 20 m	royal blue	1 piece	14826
8 cm x 20 m	royal black	1 piece	14827

Carton units: 12 pieces

size	colour	pcs./pack. unit	Item No.
10 cm x 20 m	red	1 piece	14831
10 cm x 20 m	dark blue	1 piece	14832
10 cm x 20 m	green	1 piece	14833
10 cm x 20 m	purple	1 piece	14834
10 cm x 20 m	yellow	1 piece	14835
10 cm x 20 m	royal blue	1 piece	14836
10 cm x 20 m	royal black	1 piece	14837

Carton units: 10 pieces

ADVANTAGES

- adheres only to itself and not to the skin, fur or wound
- non-slip & breathable

APPLICATION AREAS

- for fixing and relieving stress on often moved body parts



ROGGHAFT Color

Cohesive, elastic fixation bandage

This self-adhesive, non-slip, air-permeable and skin-friendly, finely woven fixation bandage is gas and steam sterilizable. It can be stretched to about 90 %. Stretched, it is about 20 m long. Particularly suitable as a supportive fixation bandage for the post-treatment of distortions.

Ingredients: 39 % cotton, 61 % polyamide, with latex coating

size	colour	pcs./pack. unit	Item No.
6 cm x 20 m	red	1 piece	14601
6 cm x 20 m	green	1 piece	14602
6 cm x 20 m	purple	1 piece	14603
6 cm x 20 m	yellow	1 piece	14604
6 cm x 20 m	royal blue	1 piece	14605
6 cm x 20 m	royal black	1 piece	14606

Carton units: 16 pieces

size	colour	pcs./pack. unit	Item No.
8 cm x 20 m	red	1 piece	14611
8 cm x 20 m	green	1 piece	14612
8 cm x 20 m	purple	1 piece	14613
8 cm x 20 m	yellow	1 piece	14614
8 cm x 20 m	royal blue	1 piece	14615
8 cm x 20 m	royal black	1 piece	14616

Carton units: 12 pieces

size	colour	pcs./pack. unit	Item No.
10 cm x 20 m	red	1 piece	14621
10 cm x 20 m	green	1 piece	14622
10 cm x 20 m	purple	1 piece	14623
10 cm x 20 m	yellow	1 piece	14624
10 cm x 20 m	royal blue	1 piece	14625
10 cm x 20 m	royal black	1 piece	14626


Carton units: 10 pieces

APPLICATION AREAS

- for fixing and bandaging on often moved body parts

ADVANTAGES

- adheres only to itself and not to the skin, fur or wound
- non-slip & breathable







benartis F Dainamic

Elastic, self-adhesive fixation bandage

The benartis F Dainamic is a self-adhesive, cohesive elastic fixation bandage for optimal compression pressure without constricting. It is sweat and water repellent and at the same time highly breathable. Transversally, the bandage can be torn off quickly and easily by hand, but has a high tear resistance in the longitudinal direction. Despite its compact material thickness, it is very light, extremely stretchable (unstretched 2m) and has optimal adhesive strength.

APPLICATION AREAS

- for fixing wound dressings and padding material
- for applying compression bandages
- for support during training and competitions

ADVANTAGES

- does not stick to the fur and skin despite its high adhesive strength
- water-repellent and breathable at the same time
- extremely light with compact material thickness and high adhesive strength
- tears off easily by hand and very easy to roll up

size	colour	pcs./pack. unit	Item No.
5 cm x 4.5 m	black	1 piece	120
5 cm x 4.5 m	red	1 piece	121
5 cm x 4.5 m	blue	1 piece	122
5 cm x 4.5 m	green	1 piece	123
5 cm x 4.5 m	yellow	1 piece	124
5 cm x 4.5 m	white	1 piece	125
5 cm x 4.5 m	pink	1 piece	126
5 cm x 4.5 m	light Blue	1 piece	127
5 cm x 4.5 m	light green	1 piece	128
5 cm x 4.5 m	purple	1 piece	129
5 cm x 4.5 m	orange	1 piece	205106

Carton units: 36 pieces, 432 pieces

size	colour	pcs./pack. unit	Item No.
7.5 cm x 4.5 m	black	1 piece	130
7.5 cm x 4.5 m	red	1 piece	131
7.5 cm x 4.5 m	blue	1 piece	205107
7.5 cm x 4.5 m	green	1 piece	133
7.5 cm x 4.5 m	yellow	1 piece	134
7.5 cm x 4.5 m	white	1 piece	135
7.5 cm x 4.5 m	pink	1 piece	136
7.5 cm x 4.5 m	light Blue	1 piece	137
7.5 cm x 4.5 m	light green	1 piece	205108
7.5 cm x 4.5 m	purple	1 piece	139
7.5 cm x 4.5 m	orange	1 piece	205109

Carton units: 24 pieces, 288 pieces, 2,880 pieces

size	colour	pcs./pack. unit	Item No.
10 cm x 4.5 m	black	1 piece	140
10 cm x 4.5 m	red	1 piece	141
10 cm x 4.5 m	blue	1 piece	142
10 cm x 4.5 m	green	1 piece	143
10 cm x 4.5 m	yellow	1 piece	144
10 cm x 4.5 m	white	1 piece	145
10 cm x 4.5 m	pink	1 piece	146
10 cm x 4.5 m	light Blue	1 piece	147
10 cm x 4.5 m	light green	1 piece	148
10 cm x 4.5 m	purple	1 piece	149
10 cm x 4.5 m	orange	1 piece	205110

Carton units: 18 pieces, 216 pieces, 2,160 pieces

size	colour	pcs./pack. unit	Item No.
15 cm x 4.5 m	blue	1 piece	152

Carton units: 12 pieces, 144 pieces, 1,440 pieces







## Valueline Flex Bandage

Colored, elastic, self-adhesive fixing bandage

The Valueline Flex Bandage is a self-adhesive, highly flexible bandage made of cotton and polyamide with a latex coating. It is non-slip, breathable and can be put on perfectly. Available in four colors. It tears easily off by hand.

Ingredients: cotton, polyamide, latex

size	colour	pcs./pack. unit	Item No.
10 cm x 4.5 m	red	10 pcs.	203598
10 cm x 4.5 m	dark green	10 pcs.	203599
10 cm x 4.5 m	blue	10 pcs.	203600
10 cm x 4.5 m	black	10 pcs.	203601

Carton units: 200 pieces

### APPLICATION AREAS

- for fixing wound dressings and padding material
- for applying compression bandages – especially on moving parts of the body

### ADVANTAGES

- available in four different colors
- tears off easily by hand
- does not stick to skin or fur
- enables painless dressing changes
- quality product from the market leader



**Supplementary product:**  
EQUIMOLL® ORIGINAL



## EQUIZINK

Creamy, moist zinc paste bandage / Unna's paste

EQUIZINK is a creamy, moist as well as elastic bandage consisting of an elastic cotton gauze bandage impregnated with a creamy matrix of water, zinc oxide and hydrocolloid. Skin-friendly and quick-drying. Ready-to-use, individually and airtight packed in a special aluminum foil bag.

Ingredients: cotton, zinc oxide, hydrocolloid, water

size	pcs./pack. unit	Item No.
10 cm x 10 m	1 pcs.	203612
10 cm x 10 m	10 pcs.	204261

Carton units: 20 pieces

### APPLICATION AREAS

- for the relief of pain and for the semi-rigid relief of injured and stressed extremities of horses with inflammatory diseases, traumatic injuries and swellings
- also suitable for effective first aid!

### ADVANTAGES

- has a cooling and analgesic effect at the same time
- the high elasticity enables the bandage to be applied evenly, lengthways and crossways, without wrinkles
- after drying, the semi-rigid zinc paste bandage calms and relieves the affected part of the body



For veterinarians only

**For different areas such as:**

- eyes
- agility
- joints
- skin and more.

A selection of our most popular products can be found starting with page 72.

We would be happy to advise you personally. Give us a call!  
+49 (0) 8133/9320-10





EQUI-SOFT

Soft padded wadding

Particularly soft cotton pad made from a cotton-viscose mixture without any chemical binding agents. Air-permeable, highly absorbent, even layer thickness, no lacings and tears off easily by hand. It can be applied quickly and easily to all parts of the body without wrinkles and is very easy to unwind.

Ingredients: cotton, viscose

size	colour	pcs./pack. unit	Item No.
20 cm x 6 m	white	1 roll	203637

Carton units: 30 rolls

size	colour	pcs./pack. unit	Item No.
38 cm x 6 m	white	1 roll	203704

Carton units: 8 rolls

APPLICATION AREAS

- for padding under various bandages;
- especially for hoof and claw bandages
- suitable for direct application on uninjured, closed skin

ADVANTAGES

- thanks to the special cotton-viscose mixture, it is particularly soft and absorbent with uniform thickness
- tears off easily by hand



benartis cotton wool Dainamic 50/50

Dressing wadding

Bandage wadding made of 50 % cotton and 50 % viscose. It is particularly soft and skin-friendly and, therefore, suitable for direct placing on closed skin. Can be applied quickly and easily on all parts of the body. It tears off by hand and the individual layers adhere well to one another and do not slip.

Ingredients: 50 % cotton, 50 % viscose

size	colour	pcs./pack. unit	Item No.
40 cm x 3 m	white	1 piece	205111

Carton units: 20 pieces, 400 pieces

APPLICATION AREAS

- for padding under various bandages; especially for hoof and claw bandages
- suitable for direct application on uninjured, closed skin



ROGG synthetic cotton bandage

Particularly soft padding bandage

Very soft, 3 mm thick cotton bandage is made of 100 % binder-free polyester with a high level of air and liquid permeability and a wrinkle-free fit.

Ingredients: 100 % binder-free polyester

size	colour	pcs./pack. unit	Item No.
6 cm x 3 m	white	1 piece	1320

Carton units: 108 pieces

size	colour	pcs./pack. unit	Item No.
8 cm x 3 m	white	1 piece	1324

Carton units: 72 pieces

size	colour	pcs./pack. unit	Item No.
10 cm x 3 m	white	1 piece	1321

Carton units: 72 pieces

size	colour	pcs./pack. unit	Item No.
12 cm x 3 m	white	1 piece	1325

Carton units: 36 pieces

size	colour	pcs./pack. unit	Item No.
15 cm x 3 m	white	1 piece	1322

Carton units: 36 pieces

APPLICATION AREAS

- for padding hard edges of splint bandages
- for padding under compression bandages

ADVANTAGES

- high air and moisture permeability



benartis cotton wool Xactic

Dressing wadding

Padding bandage made of 100 % synthetic material – non-crinkled and non-absorbent polyester fibers. It is particularly soft and sturdy. Can be unwound very easily and, due to its skin-friendly properties, is suitable for direct application on closed skin. It can be applied quickly and easily to keep its shape. Due to its fiber structure, the individual bandage sections adhere well to one another and do not slip.

size	colour	pcs./pack. unit	Item No.
10 cm x 3 m	white	1 piece	196

Carton units: 120 pieces, 2,880 pieces

APPLICATION AREAS

- for padding under various bandages; especially for hoof and claw bandages
- suitable for direct application on uninjured, closed skin





ROGG foam bandage

Thin underlay bandage, 4 mm

Ideal padding and underlay material for support, compression and rigid bandages.

Ingredients: 100 % polyurethane foam

size	colour	pcs./pack. unit	Item No.
10 cm x 2 m	white	1 piece	202648
15 cm x 2 m	white	1 piece	202649

APPLICATION AREAS

- for padding under various bandages
- suitable for direct application on uninjured, closed skin

ADVANTAGES

- extra thin
- high cushioning effect



Supplementary product:  
EQUIMOLL® ORIGINAL



snögg AnimalPOLSTER

Cross-elastic, self-adhesive  
foam bandage – 5 mm

The cross-elastic, soft foam bandage, made of polyurethane foam, with an adhesive side, is with its thickness of 5 mm quite robust. It is porous and, therefore, breathable, absorbent and adheres very well to fur / hair and skin, but it can still be removed effortlessly and painlessly. No direct contact to wounds!

size	colour	pcs./pack. unit	Item No.
14 cm x 2 m	white	1 piece	203743

Carton units: 30 pieces

size	colour	pcs./pack. unit	Item No.
9 cm x 2 m	white	1 piece	204968

Carton units: 48 pieces

APPLICATION AREAS

- for padding and fixing wound dressings and bandages – especially on difficult parts of the body

ADVANTAGES

- latex free
- good adhesive properties
- high protection thanks to robust padding

Dental tartar in dogs and cats is one frequent diagnosis – in your practice too?

RECOMMEND DERMAVET LAVAPIROX® VET ANTI-PLAQUE C TO YOUR PATIENT OWNERS FOR DAILY ORAL HYGIENE!

DERMAVET LAVAPIROX® Vet Anti-Plaque C – developed and manufactured by veterinarians in Germany – prevents gum inflammation, bad breath and tooth loss and spares the animal chronic pain.

**Mode of action:** DERMAVET LAVAPIROX® Vet Anti-Plaque Anti-Plaque C contains chlorhexidine to combat bacteria in the oral cavity. Thus, with regular use, there is no biofilm on the teeth and no bad breath. The contained specially adapted micro-abrasive particles guarantee a good cleaning effect without damaging the enamel, even with short cleaning times.



ORDER NOW AND BENEFIT FROM OUR FAVORABLE RESALE CONDITIONS!

DERMAVET LAVAPIROX®  
Anti-plaque C

Toothpaste for dogs and cats

size	pcs./pack. unit	Item No.
50 ml	1 piece	205084

Carton units: 12 pieces

APPLICATION AREAS

- for daily teeth cleaning

ADVANTAGES

- effective against plaque bacteria
- prevents inflammation, thanks to chlorhexidine





EQUI-PLAST

Elastic, adhesive plaster bandage / adhesive bandage

EQUI-PLAST is a white, air-permeable, waterproof and extremely resistant adhesive bandage, made of 100 % cotton, coated on one side with zinc oxide rubber adhesive, without paper. The extra strong adhesive bandage is 75 % lengthways elastic (stretched approx. 4.5m long), does not slip off even during movement and can be removed very easily without pulling on the fur.

Ingredients: cotton, zinc oxide rubber glue

size	colour	pcs./pack. unit	Item No.
10 cm x 4.5 m	white	1 roll	203613

Carton units: 50 pieces

ADVANTAGES

- extremely strong adhesive behavior despite 75 per-cent longitudinal elasticity

APPLICATION AREAS

- first aid for minor injuries and abrasions
- fixation of wound dressings for larger wounds
- as a mobilization compression bandage for muscle and ligament injuries
- in combination with tape pulls
- as an elastic substitute tape



benartis Plaster Bandage

Stretchable bandage

This adhesive and elastic bandage consists of a woven cotton backing and does not contain any natural latex. The elasticity and length guarantee a perfect fit even during movement. The bandage is permeable to air and allows the skin to breathe. Due to its special structure, the bandage can be rolled up gradually, and thus more easily and precisely applicable to joints or muscles. Ideal for bandages where a secure hold is required. This is especially the case with stronger compression bandages, sling bandages and particularly in equestrian sports.

size	colour	pcs./pack. unit	Item No.
10 cm x 4.5 m (stretched)	white	1 piece	1150

Carton units: 12 pieces, 144 pieces, 1,440 pieces

APPLICATION AREAS

- For bandages with particularly strong hold such as:
- stronger compression bandages
  - sling dressings
  - ideal for equestrian sports



ROGG Stretch

Self-adhesive, cut-to-size fixation bandage

White, cross-elastic polyester fleece with polyacrylate adhesive. Micro-perforated with high adhesion on shaved hair / fur.

Ingredients: polyester, polyacrylate adhesive, fleece

size	colour	pcs./pack. unit	Item No.
10 cm x 10 m	white	1 piece	1661
15 cm x 10 m	white	1 piece	1662
20 cm x 10 m	white	1 piece	1663

Carton units: 36 pieces

APPLICATION AREAS

- as a secondary bandage or secondary plaster for the fixation of wound dressings

ADVANTAGES

- easy to cut
- does not fray
- high adhesion on shaved hair / fur

**Supplementary product:**  
Zemuko® cellulose non-woven compress







Miro Plast

Sterile, elastic first-aid wound dressing

Miro Plast Steril is a breathable, hypoallergenic and sterile wound dressing for the treatment of injuries, cuts and tears. The wound pad made of polyester mesh is particularly absorbent and thus, adhesion to the wound is prevented. The very skin-friendly, hypoallergenic polyacrylate adhesive can be gently removed from the shaved skin without leaving any adhesive residue.

Ingredients: terylene, viscose, polyacrylate adhesive, polyester

size	colour	pcs./pack. unit	Item No.
5 cm x 7.2 cm	white	100 pcs.	202001
10 cm x 8 cm	white	50 pcs.	202002
15 cm x 10 cm	white	50 pcs.	202003
20 cm x 10 cm	white	50 pcs.	202004

APPLICATION AREAS

- for the treatment of general wounds – but especially cuts and tears.

ADVANTAGES

- no sticking to the wound and high absorbency at the same time



ROGG Silk

White universal plaster with a zigzag edge

This air-permeable universal plaster, made of taffeta fabric, can be torn lengthways and crossways. It can be easily written on. The special zinc oxide polyacrylate adhesive does not contain any solvents and can be removed again without leaving any residue.

Ingredients: 100 % cotton and cellulose fibers, polyacrylate glue

size	colour	pcs./pack. unit	Item No.
2.5 cm x 9 m	white	12 pieces	1631

APPLICATION AREAS

- for secure bandage fixation

ADVANTAGES

- tears off by hand
- solvent free

OUR RECOMMENDATION FOR TYPICAL EQUINE SKIN DISEASES

FOR SWEET ITCH



DERMAVET LAVAPIROX® foam against fungi, bacteria and indeterminate itching.

Equital Immune Active to support the immune system, wound healing and well-being.

FOR SCURF



DERMAVET LAVAPIROX® shampoo and foam for weekly and daily fur cleaning. Accompanying therapy against the often difficult to determine triggers of scurf. Dissolves dandruff and penetrates crusts. Particularly gentle due to a pH value of 7.




IN ADDITION: FOR RINSING AND DISINFECTION

DERMAVET LAVAPIROX® Fluid helps you treat fungi and bacteria, especially in hard-to-reach places, and disinfects without burning.

You can find out more about our care products on pages 72 – 79.





## EQUICOMP

### Elastic tubular bandage

Lengthways and crosswise elastic tubular bandage made of 100 % cotton, seamlessly knitted. It can be cut at any point, is run-proof and has a good anatomical, crease-free fit.

Ingredients: cotton

size	pcs./pack. unit	Item No.
12 cm x 5 m	1 roll loose	203640


Carton units: 60 rolls loose

### APPLICATION AREAS

- as a pressure-free bandage
- as an under- and overcoat for zinc paste, cast and plaster casts as well as for splint and compression bandages
- for fixing wound dressings and padding material

### ADVANTAGES

- no pressure points due to the seamless processing
- easy handling due to high elasticity
- well absorbent and air permeable at the same time due to the 100 % cotton



## ROGG terry cloth tubes

### Universal padding and underlay material

This terry cloth tube made of soft natural fibers is 10 m long, knitted seamlessly and can be applied without creases. It can be cut at any point, is run-proof and, thanks to its precise fit, is extremely comfortable to wear.

Ingredients: 98 % cotton, 2 % elastane

size	colour	pcs./pack. unit	Item No.
L - 40-75 cm	white	1 piece	88105
L - 40-75 cm	pink	1 piece	88115
L - 40-75 cm	blue	1 piece	88125
L - 40-75 cm	green	1 piece	881350

XL - 50-95 cm	white	1 piece	88106
XL - 50-95 cm	pink	1 piece	88116
XL - 50-95 cm	blue	1 piece	88126
XL - 50-95 cm	green	1 piece	881360

### APPLICATION AREAS


- as a pressure-free bandage
- as an under- and overcoat for zinc paste, cast and plaster casts as well as splint and compression bandages
- for fixing wound dressings and padding material

### ADVANTAGES

- no pressure points due to the seamless processing
- easy handling due to high elasticity



The basis for your treatment regimen according to the principles of modern wound care: not only protect wounds with a bandage, but also use modern wound antiseptics and special wound dressings to greatly accelerate wound healing in all stages, and to combat infections in a targeted and gentle manner. Objective: to create a warm, moist, infection-free wound environment in order to provide maximum support for the cells in their healing phase.







KerraSol™ wound irrigation solution

Ecological wound cleansing solution  
– water and salt based

Basic product for cleaning, moistening and infection prophylaxis in the wound and for removing encrusted wound coatings, bandages and coverings. KerraSol™ wound irrigation solution contains only natural, endogenous ingredients, is ready to use and has a particularly long shelf life (at least 2 years unopened, after opening at least 12 weeks). It offers painless application, an effective reaction in the wound, is effective against Pseudomonas and eliminates wound odors.

Ingredients: water, sodium hypochlorite <0.06 % (electrochemically activated mineral salt solution)

size	pcs./pack. unit	Item No.
250 ml	1 piece	203539
1000 ml	1 piece	205087

Carton units: 6 pieces

APPLICATION AREAS

- for cleaning & humidification
- for antimicrobial decontamination and infection prophylaxis of wounds
- for the painless removal of wound dressings

ADVANTAGES

- no surfactants, no alcohol
- no restriction in the duration of use
- deposits are removed
- wound odors are eliminated quickly and reliably
- non-irritating
- very well tolerated
- painless



KerraSol™ wound spray gel

Ecological hydrogel in a spray form  
– water and salt based

The KerraSol™ wound spray gel is easy to apply and remains in the wound as a hydrogel. It is used for long-term moistening and moisture retention of wounds, forms an ideal wound healing environment and, thus, supports the natural wound healing process. The preservative (sodium hypochlorite) contained has a very broad spectrum against bacteria (including Pseudomonas aeruginosa, MRSA / MRE), viruses and fungi, which enables non-specific decontamination of wounds.

Ingredients: water, sodium hypochlorite <0.06 % (electrochemically activated mineral salt solution), lithium-magnesium-sodium-silicate

size	pcs./pack. unit	Item No.
75 ml	1 piece	203540

APPLICATION AREAS

- for cleaning and long-term humidification
- for antimicrobial decontamination and Infection prophylaxis of wounds
- for the painless removal of wound dressings

ADVANTAGES

- no surfactants, no alcohol
- no restriction in the duration of use
- biofilm and fibrinoid deposits are removed
- wound odors are eliminated quickly and reliably
- non-irritating
- very well tolerated
- painless



ROGG Bio-Kollagen (organic collagen)

Sterile, fully absorbable collagen  
wound dressing

This fleece-like wound filler, that is completely absorbable by the body, supports hemostasis, lowers the pH value in the wound environment and accelerates wound closure. It can be easily modeled on the wound size, does not swell in the wound, does not stick to the wound base and, thus, avoids pain. Also suitable for operative tissue defects, wound cavities and burns.

Ingredients: 100 % collagen, from the skin of certified Demeter pigs

size	colour	pcs./pack. unit	Item No.
2.2 cm x 2.7 cm	white	5 pieces	200998
2.2 cm x 2.7 cm	white	10 pieces	200808
4.5 cm x 5.5 cm	white	5 pieces	200999
4.5 cm x 5.5 cm	white	10 pieces	200809
4.5 cm x 11 cm	white	5 pieces	200997
4.5 cm x 11 cm	white	10 pieces	200807
9 cm x 11 cm	white	5 pieces	201000
9 cm x 11 cm	white	10 pieces	200801

APPLICATION AREAS

- for wound filling in deep wounds
- suitable for acute and deep wounds in all wound healing phases

ADVANTAGES

- 100 % absorbed by the body



ROGG Vulcosan AMC

Sterile, polyhexanide, moist compress

This MRSA-effective, PHMB\*-impregnated wound dressing is used for effective decontamination and infection prophylaxis against microorganisms. The contained anti-microbial complex (AMC) dissolves the biofilm, cleans the wound and creates an optimal wound environment for the formation of new granulation tissue – without maceration of the wound edge.

Ingredients: viscose, \*0.1 % polyhexamethylene biguanide

size	colour	pcs./pack. unit	Item No.
10 cm x 10 cm	white	10 pieces	201467
10 cm x 10 cm	white	50 pieces	203351

APPLICATION AREAS

- supports the healing of acute and chronic wounds with an increased risk of infection
- especially suitable for superficial, deep and / or post-operative wounds and burns up to 2nd degree
- for packing
- for use under compression

ADVANTAGES

- prevents biofilm
- disinfects and decontaminates wounds at risk of infection





curea clean Veterinary

Superabsorbent secondary wound dressing

Wound dressing with hydroactive surface and special SuperCore® absorbent core for direct application to wounds with little to very high levels of exudation. The cellulose fleece on the wound side serves to remove secretions. The special absorbent core takes up the secretion along with germs and cell debris. The absorption capacity is 75 g / 100 cm². Without latex, adhesives or binders.

Ingredients: polypropylene, cellulose, sodium polyacrylate

size	colour	pcs./pack. unit	Item No.
10 cm x 10 cm	light blue	100 pieces	203719
10 cm x 20 cm	light blue	50 pieces	203722
15 cm x 15 cm	light blue	50 pieces	203720
20 cm x 20 cm	light blue	25 pieces	203721
20 cm x 30 cm	light blue	25 pieces	203723

APPLICATION AREAS

- wound dressing for low to heavily exuding wounds

ADVANTAGES

- adhesive- and binder-free, cellulose-based absorbent core
- high absorption capacity and speed even under pressure
- secure retention
- reliable binding of germs
- uniform liquid distribution in the absorbent core
- high flexibility
- strong odor binding
- minimized risk of wound maceration



Revamil® wound dressing

Honey wound compress / pad

This wound compress is impregnated with pure honey produced under controlled conditions. The Revamil® honey pad is of reproducible quality and guarantees no residues of pesticides. The honey contains a high content of enzymes which kill bacteria effectively and quickly. With its high sugar concentration, the honey combines with the wound fluid and forms a reinforcing, antibacterial protective layer of the compress. Thanks to the hydrophilic properties of Revamil® honey and its low pH value, a moist and favorable wound environment is created for faster, natural wound healing.

Ingredients: medical honey, pulp, sodium polyacrylate, non-woven polypropylene, polyethylene

size	colour	pcs./pack. unit	Item No.
5 cm x 5 cm	light yellow	10 pieces	202783
8 cm x 8 cm	light yellow	7 pieces	202784
10 cm x 20 cm	light yellow	5 pieces	202785

APPLICATION AREAS

- decontaminates and protects wounds from infections with the help of antibacterial properties
- particularly suitable for 1st and 2nd degree burns
- for the treatment of chronic, infected, acute or surgical wounds

ADVANTAGES

- forms a reinforcing, antibacterial protective layer between the wound and the compress
- optimizes the pH value of the wound environment thanks to its hydrophilic properties



Revamil® gel

Hydrophilic honey wound gel

This wound gel contains pure honey produced under controlled conditions, which has a reproducible quality and guarantees no residual levels of pesticides. With its high content of enzymes, Revamil® honey kills bacteria effectively and quickly. Due to its high sugar concentration, the honey combines with the wound fluid and thus forms a reinforcing, antibacterial protective layer. Thanks to the hydrophilic properties of Revamil® honey and its low pH value, a moist and favorable wound environment is created for faster, natural wound healing. The single-dose syringe brings the honey wound gel well into deep wounds and under the wound edges.

Ingredients: medical honey

size	colour	pcs./pack. unit	Item No.
2 g single dose	light yellow	12 pieces	202786
18 g tube	light yellow	1 piece	202787

APPLICATION AREAS

- decontaminates and protects wounds from infections thanks to the antibacterial effect
- particularly suitable for the treatment of chronic wounds, infected, acute, surgical wounds or first and second degree burns

ADVANTAGES

- forms a reinforcing, antibacterial protective layer
- optimizes the pH value of the wound environment thanks to its hydrophilic properties



Revamil® balm

Hydrophilic ointment with honey

The easy-to-apply ointment contains pure honey produced under controlled conditions, which has a reproducible quality and guarantees no residual levels of pesticides. With its high content of enzymes, Revamil® honey kills bacteria effectively and quickly. Due to its high sugar concentration, the honey mixes with the wound fluid and thus forms a reinforcing, antibacterial protective layer. Thanks to the hydrophilic properties of Revamil® honey and its low pH value, a moist and favorable wound environment is created for faster, natural wound healing. Without preservative substances.

Ingredients: peanut oil, beeswax, oleic acid glycerol ester, medicinal honey, water

size	colour	pcs./pack. unit	Item No.
15 g tube	light yellow	1 piece	202788

APPLICATION AREAS

- for the treatment of superficial wounds

ADVANTAGES

- forms a reinforcing, antibacterial protective layer
- optimizes the pH value of the wound environment thanks to its hydrophilic properties





### NAWAlution Gel

#### Odorless hydro wound gel

The alcohol-free, hydroactive wound gel, effective against MRSA, cleans and moisturizes wounds gently and, at the same time, decontaminates dry and stubborn wound coverings. By optimizing the pH value and, thus, the wound environment, the chronic, stagnant healing phase can be overcome. NAWAlution Gel is suitable for repeated and long-term use. Can also be used on mucous membranes.

Ingredients: wificata, hydroxypropyl starch phosphate, cocamidopropyl betaine, trace elements (zinc, iron), polyaminopropyl biguanide (PHMB)

size	pcs./pack. unit	Item No.
5 ml	10 pieces	40060

#### APPLICATION AREAS

- for gentle cleaning and effective decontamination of dry, stubbornly covered wounds
- promotes healing
- for immediate cooling of 1st and 2nd degree burns
- for moistening bandages and wound dressings

#### ADVANTAGES

- moistens, disinfects and decontaminates necrotic wounds
- optimizes the pH value of the wound environment



### ROGG wound & burn ointment

#### Cooling, versatile ointment

This ointment with a cooling effect relieves pain, protects the wound bed from drying out, prevents germ penetration and creates an optimal climate for wound healing. It promotes the healing of injured skin, which results in accelerated regeneration of tissue and significantly reduced scarring. The ointment does not contain any skin-irritating substances and can also be used in the mouth and mucous membrane area.

Ingredients: aqua purificata, octyldodecanol, hydrogenated palm glycerides, PEG-200 hydrogenated glyceryl palmate, cetyl alcohol, methylparaben, propylparaben, sorbic acid, dehydroacetic acid, trace elements (zinc, iron)

size	pcs./pack. unit	Item No.
100 g	1 piece	36021

Carton units: 31 pieces

#### APPLICATION AREAS

- for setting an optimal physical environment for the healing of injured skin
- well suited for superficial wounds, cuts, abrasions, tears as well as bruises, contusions, and 1st and 2nd degree burns



### snögg AnimalVET

#### Transversally elastic, self-adhesive wound foam dressing – 5 mm

Cross-elastic, self-adhesive padding made of very robust, flexible, absorbent foam. Does not slip, does not absorb. Can be used as a wound pad. Stops bleeding. Stays in place even in water.

size	colour	pcs./pack. unit	Item No.
8 cm x 5 m	white	1 piece	204966
8 cm x 1 m	white	1 piece	205103
14 cm x 5 m	white	1 piece	205104


Carton units: 20 pieces

#### APPLICATION AREAS

- for dressing, padding and fixation of injuries and open wounds
- for initial dressing

#### ADVANTAGES

- adheres only to itself and not to skin, fur or wound



### snögg AnimalSOFT

#### Transversally elastic, self-adhesive wound foam dressing – 1.9 mm

Cross-elastic, self-adhesive padding made of extra thin, flexible, absorbent foam. Does not slip, does not absorb. Can be used as a wound pad. Stops bleeding. Stays in place even in water.

size	colour	pcs./pack. unit	Item No.
6 cm x 4.5 m	white	1 piece	205105

Carton units: 20 pieces

#### APPLICATION AREAS

- for dressing, padding and fixation of injuries and open wounds
- for initial dressing

#### ADVANTAGES

- extra thin and therefore extra flexible
- does not restrict movement
- adheres only to itself and not to skin, fur or wound





## ROGG washing lotion

## Nourishing washing lotion for stressed skin

This washing lotion with skin-friendly pH value and high-quality ingredients gently cleanses the skin and cares for it at the same time. It preserves the physiological skin environment and prevents the skin from drying out. It is very economical in use.

☐ Ingredients: aqua, sodium laureth sulfate, sodium chloride, PEG-7 glyceryl, cocoate, cocamide DEA, sodium cocosulfate, parfum, 2-bromo-2-nitropropane-1,3-diol, phenoxyethanol, citric acid, CI 42051

size	pcs./pack. unit	colour	Item No.
500 ml dispenser bottle	1 piece	blue	9712

 Carton units: 20 pieces

size	pcs./pack. unit	colour	Item No.
1000ml dispenser bottle	1 piece	blue	9711

 Carton units: 12 pieces

size	pcs./pack. unit	colour	Item No.
5l canister	1 piece	blue	9710

## APPLICATION AREAS

- for hand washing

## ADVANTAGES

- mild cleaning and care with skin-friendly PH value



Highly effective hand disinfectant

The skin/hand disinfectant has been filtered free of spores, is bactericidal (incl. TbB) and fungicidal / levurocidal. Do not use on mucous membranes or in the eye area.

Ingredients: 2-propanol, purified water

size	pcs./pack. unit	colour	Item No.
500ml dispenser bottle	1 piece	white	9600

 Carton units: 20 pieces

size	pcs./pack. unit	colour	Item No.
1l dispenser bottle	1 piece	white	9601

 Carton units: 12 pieces

size	pcs./pack. unit	colour	Item No.
5l canister	1 piece	white	9602

## APPLICATION AREAS

- for hygienic and surgical hand disinfection



## ROGG hand disinfectant gel

Fragrance-free, with a moisturizing effect

The formulation of this hand disinfectant as a gel enables the use of a moisturizing substance, so that the hands are not damaged despite frequent disinfection. It is bactericidal (including MRSA and tuberculosis), fungicidal / levurocidal (*Candida albicans*) and virucidal (including noroviruses and corona viruses).

 Ingredients: ethanol, water, gelling agent, moisturizing agent

size	pcs./pack. unit	colour	Item No.
150 ml gown dispenser bottle	1 piece	white	9605
500 ml dispenser bottle	1 piece	white	9603

 Carton units: 20 pieces

size	pcs./pack. unit	colour	Item No.
1l dispenser bottle	1 piece	white	9604

 Carton units: 12 pieces

## APPLICATION AREAS


- for hygienic hand disinfection and care
- ROGG hand disinfectant gel

## ADVANTAGES

- moisturizing effect thanks to gel formulation

OUR SKIN & HAND DISINFECTION IS  
**EFFECTIVE AGAINST:**

- Hepatitis B Virus (HBV)
- Hepatitis C virus (HCV)
- Herpes viruses
- Bovine Viral Diarrhea Virus (BVDV)
- Human Immunodeficiency Virus (HIV)
- Influenza A, Influenza B
- Methicillin resistant
- Staphylococcus aureus (MRSA)
- Noroviruses (MNV)
- Vaccinia viruses (VACV / VV)
- Salmonella
- Rotaviruses
- Corona viruses







Isopropyl alcohol 70 %

Disinfectant & cleaning agent

Isopropyl alcohol 70 % is a disinfectant for hygienic and surgical hand disinfection, skin disinfection before simple injections and punctures of peripheral vessels, skin disinfection prior to operations and before punctures of joints, disinfection of skin rich in sebum and for the preparation of cooling compresses.

Ingredients: 2-propanol, purified water

size	pcs./pack. unit	colour	Item No.
1l	1 bottle	colorless	1691
5l	1 canister	colorless	16901

! Always use isopropyl alcohol 70 % exactly according to the package insert.



Isopropyl alcohol > 99 %

Disinfectant & cleaning agent

Isopropyl alcohol > 99 % is highly concentrated alcohol and can be used as a substitute for ethanol. It is suitable for degreasing, disinfecting and cleaning many material surfaces. The product is readily soluble in water, evaporates quickly and is residue-free. It can also be used for cleaning and as a solvent for various organic substances.

Ingredients: 2-propanol

size	pcs./pack. unit	colour	Item No.
5l	1 canister	colorless	204914
10l	1 canister	colorless	204909

! Always use isopropyl alcohol 99.9 % exactly according to the package insert.



kodan® tincture forte, colorless

Alcoholic antiseptic for application on the skin

For preoperative skin antisepsis and postoperative suturing. Before catheterizations, injections, blood sampling, punctures, excisions, cannulations and biopsies. To support general hygienic measures in the prevention of fungal skin infections and for hygienic hand disinfection. The product is not suitable for mucous membranes and larger wounds.

Ingredients: propan-2-o, propan-1-ol, biphenyl-2-ol, hydrogen peroxide solution 30 % (Ph. Eur.), purified water

size	colour	pcs./pack. unit	Item No.
250 ml	colorless	1 bottle	99001

Carton units: 10 bottles

size	colour	pcs./pack. unit	Item No.
1l	colorless	1 bottle	204671
2l	colorless	1 bottle	9899

APPLICATION AREAS

- for skin disinfection before surgical interventions, catheterization, blood and liquor sampling, injections, punctures, excisions, cannulations, biopsies, etc.
- for suturing
- to support general hygienic measures
- as part of the prevention of fungal skin diseases

+ Supplementary product: ROGG gauze swabs



kodan® tincture forte, colored

size	colour	pcs./pack. unit	Item No.
250 ml	brown	1 bottle	99004
1l	brown	1 bottle	204670
2l	brown	1 bottle	9900

Ingredients: propan-2-o, propan-1-ol, biphenyl-2-ol, hydrogen peroxide solution 30 % (Ph. Eur.), purified water, quinoline yellow, yellow orange S, brilliant black BN





## ROGG disinfectant wipes AF Q\*

For fast, effective cleaning and disinfection of objects

These aldehyde- and phenol-free disinfectant wipes have been dermatologically tested, are well tolerated by the skin and dry without leaving residues. They are bactericidal (including TbB / MRSA), fungicidal / yeasticidal and virus inactivating. Suitable for the disinfection of surfaces, loungers and parts made of plastic such as ultrasound heads or plexiglass.

Ingredients: didecyldimethylammonium chloride, quaternary ammonium compounds, benzyl C12-16 alkyl dimethyl chloride

size	colour	pcs./pack.	unit	Item No.
Dispenser box	white	100 wipes		9743

Carton units: 12 dispenser cans

size	colour	pcs./pack.	unit	Item No.
Refill pack	white	100 wipes		9744

### APPLICATION AREAS

- for quick, effective cleaning and disinfection of objects

\*AF = aldehyde-free



## ROGG rapid disinfection set AF plus\*\*

Ready to use and optimally dosed

The ROGG quick disinfection set AF for 3 fillings can be created by ordering: 1x ROGG system dispenser bucket, 3x 2l ROGG rapid disinfection AF plus as well as 1x 3 rolls of ROGG system non-woven wipes with 100 sheets each. The dry ROGG system fleece cloths consist of high quality, lint-free spunlace perforated fleece, knitted from a 100% synthetic fiber blend. They can be impregnated with any surface disinfectant and optimally release the active ingredient onto the surface. 1 roll consists of 100 wipes á 32 cm x 30 cm.

### ROGG system dispenser bucket (order 1x)

colour	pcs./pack.	unit	Item No.
blue	1 piece (empty)		9761

### ROGG system non-woven wipes (order 1x)

size	colour	pcs./pack.	unit	Item No.
100 sheets / roll	white	3 roles		203529

### ROGG rapid disinfection (order 3x)

size	colour	pcs./pack.	unit	Item No.
2l dispenser bottle	white	1 piece		9763

### APPLICATION AREAS

- for disinfecting surfaces

\*\*AFQ = aldehyde-free, alcohol-free



## ROGG system non-woven wipes

Refill pack for ROGG rapid disinfection set AF


The dry ROGG system non-woven wipes consist of high quality, lint-free spunlace perforated fleece, knitted from a 100% synthetic fiber blend. They can be impregnated with any surface disinfectant and optimally release the active ingredient onto the surface.

1 roll consists of 100 wipes á 32 cm x 30 cm.

size	colour	pcs./pack.	unit	Item No.
100 sheets	white	3 roles		203529
100 sheets	white	1 roll		205006

### APPLICATION AREAS

- for refilling the ROGG rapid disinfection set AF







## ROGG rapid disinfection AF plus

Ready-to-use, highly effective rapid disinfection with a short exposure time

The aldehyde- and phenol-free, dermatologically tested disinfectant is bactericidal (incl. TbB), fungicidal / yeas-ticidal and virus-inactivating. It is suitable for the rapid disinfection of medical inventory and for spraying inac-cessible areas.

Ingredients: ethanol, didecyldimethylammonium chloride

size	colour	pcs./pack. unit	Item No.
1l dispenser bottle	white	1 piece	9723
2l dispenser bottle	white	1 piece	9763
5l canister	white	1 piece	9724

### APPLICATION AREAS

For disinfecting surfaces, by wiping and spraying

### ADVANTAGES

- short exposure time
- even inaccessible areas can be reached by spraying



## ROGG rapid disinfection AF Q

Ready-to-use, highly effective rapid disinfection without alcohol

The aldehyde-, phenol- and alcohol-free disinfectant is bactericidal including MRSA, odorless and can be ap-plied as a foam. It is suitable for the rapid disinfection and cleaning of medical inventory, medical products and surfaces of all kinds, in accordance with the biocide direc-tive, especially for alcohol-sensitive surfaces and acrylic glass. For hard-to-reach areas: with foam spray head (can be applied as a foam with a neutral odor; reduction of aerosols.)

Ingredients: quaternary ammonium compounds, benzyl-C12-16-alkyldimethyl-chloride; didecyldimethy-lammonium chloride

size	colour	pcs./pack. unit	Item No.
1l dispenser bottle (without foam spray head)	white	1 piece	97275
1l dispenser bottle (including foam spray head)	white	1 piece	97271

size	colour	pcs./pack. unit	Item No.
2l bottle	white	1 piece	97272
5l canister	white	1 piece	97273

### APPLICATION AREAS

- for disinfecting surfaces, by wiping and spraying

### ADVANTAGES

- odorless
- can be applied as a foam with a foam spray head
- reduction of aerosols



## ROGG surface disinfection AF

Highly effective surface disinfectant

This aldehyde-free disinfectant is particularly user-friend-ly due to its short exposure times, high cleaning power and odorlessness. It is bactericidal (including TbB), fun-gicidal / yeas-ticidal and virus inactivating. The AF disin-fectant is used to disinfect and clean medical products, surfaces and inventory in the veterinary field, especially PVC and linoleum floors.

Ingredients: quaternary ammonium compounds, benzyl-C12-16-alkyldimethyl, chlorides

size	colour	pcs./pack. unit	Item No.
5l canister	white	1 piece	9722

### APPLICATION AREAS

- for effective disinfection and cleaning of surfaces and inventory

### ADVANTAGES

- particularly user-friendly thanks to the short exposure times
- high cleaning power
- odorless



## ROGG surface disinfection AF Q

Highly effective surface disinfectant

This aldehyde-free disinfectant is particularly user-friend-ly due to its short exposure times, high cleaning power and odorlessness. It is bactericidal (including TbB), fun-gicidal / yeas-ticidal and virus inactivating. The alcohol-free AF Q disinfectant can be used like AF, but is also suitable for disinfecting alcohol-sensitive surfaces such as acrylic glass. The contained surfactants are biodegradable accord-ing to the detergent regulation.

Ingredients: didecyldimethylammonium chloride, quaternary ammonium compounds, benzyl- C12-16-al-kyldimethyl-, chlorides.

size	colour	pcs./pack. unit	Item No.
5l canister	white	1 piece	9728

### APPLICATION AREAS

- for effective disinfection and cleaning of surfaces, inventory and acrylic glass

### ADVANTAGES

- particularly user-friendly thanks to the short exposure times
- high cleaning power
- odorless
- biodegradable





## ROGG instrument disinfection AF

Disinfection and cleaning  
with a broad range of efficacy

This aldehyde-free, material-friendly and highly effective disinfectant is bactericidal (incl. TB), fungicidal and virus-inactivating and is suitable for ultrasonic cleaning devices and waste trays.

100g of solution contain: alkylpropylenediaminguani-  
dinediacetate 78 % strength 3.2g; dialkyloxethylmethy-  
lammonium propionate 70 % strength 5.36g

size	colour	pcs./pack. unit	Item No.
5l canister	green	1 piece	9720

### APPLICATION AREAS

- for cleaning and disinfecting all types of surgical instruments, thermolabile materials and endoscopes



Supplementary product:  
Disinfection tub



## ROGG instrument disinfection VV

Particularly fast acting,  
virucidal, liquid concentrate

The formaldehyde-free formulation of this disinfectant is characterized by economical application concentrations and a broad range of efficacy with good material compatibility. It is ideal for the virucidal final disinfection of flexible and rigid endoscopes as well as medical instruments of all kinds. It is bactericidal (including MRSA) / tuberculocidal, yeastcidal and virucidal.

100g solution contain: 12g glutaral, anionic surfac-  
tants 5-15 %

size	colour	pcs./pack. unit	Item No.
2l bottle	green	1 piece	203490

### APPLICATION AREAS

- for the disinfection of endoscopes and surgical instruments of all kinds



Supplementary product:  
Disinfection tub



## Disinfection tray

Including sieve insert & lid

Filling volume: 1 liter

size	colour	pcs./pack. unit	Item No.
22.3 cm x 14.9 cm x 9.1 cm	white	1 piece	97301

Filling volume: 3 liters

size	colour	pcs./pack. unit	Item No.
20.5 cm x 16 cm x 10 cm	white	1 piece	97303

Filling volume: 4 liters

size	colour	pcs./pack. unit	Item No.
34.3 cm x 21.3 cm x 11.5 cm	transp.	1 piece	9730

Filling volume: 5 liters

size	colour	pcs./pack. unit	Item No.
39.4 cm x 26 cm x 15.6 cm	white	1 piece	97305

Filling volume: 10 liters

size	colour	pcs./pack. unit	Item No.
69.5 cm x 17 cm x 12.8 cm	white	1 piece	97310

### APPLICATION AREAS

- for device disinfection



Supplementary product:  
ROGG instrument disinfection



## Application aids

Suitable for our disinfection products

### Dosing bottle

size	colour	pcs./pack. unit	Item No.
1l dosing bottle	white	1 piece	9731

### Hand spray head for spray disinfection

size	colour	pcs./pack. unit	Item No.
for the 1l bottle	white	1 piece	9732

### Dosing pump

size	colour	pcs./pack. unit	Item No.
for the 1l & 2l bottles	white	1 piece	9739

### Canister pump

size	colour	pcs./pack. unit	Item No.
30 ml / stroke	white	1 piece	9734

### Canister tap

size	colour	pcs./pack. unit	Item No.
for the 5l canisters	white-red	1 piece	9733





Dr. Ulrich Wendlberger is a specialist veterinarian for small animals and dermatology with many years of experience in dentistry and in the treatment of typical and atypical diseases of the skin.

With DERMAVET LAVAPIROX® he has developed a product line that combats the most common symptoms and causes of itching, allergic reactions as well as fungal and bacterial infections.

35 years of professional practice and his own experience as owner of an active little terrier contributed to the development of DERMAVET LAVAPIROX®.

DERMAVET LAVAPIROX® has been designed especially for small animals such as dogs and cats. However, the DERMAVET LAVAPIROX® system has also proven effective in therapy-accompanying applications for horses, such as scurf, sweet itch and skin irritations.

## THE POWERFUL TRIO AGAINST FUNGI, BACTERIA AND ITCHING



1 2 3

DERMAVET LAVAPIROX® Shampoo, Fluid and Foam are the powerful antibacterial, antifungal team that effectively supports your treatment against the secondary causes of itching. All without negatively impacting irritated skin or wound healing. Thanks to a pH of 7, the shampoo, foam and fluid are particularly gentle on the skin of small and large animals.

### OUR APPLICATION RECOMMENDATION:

- 1 DERMAVET LAVAPIROX® Shampoo**  
Accompanying therapy for weekly coat washing in case of fungal and bacterial infestation as well as itching, dandruff and odors.
- 2 DERMAVET LAVAPIROX® Foam**  
As a daily complement for sustained bactericidal, fungicidal and keratolytic treatment. Also for waterless cleaning of dirty coats, against odors, and for perineal hygiene.
- 3 DERMAVET LAVAPIROX® Fluid**  
For spot application and rinsing, especially in hard-to-reach areas and on mucous membranes. Also for local disinfection and as an antibacterial rinse for the inner ear and foreskin.



## STRONG AND GENTLE: CLEANING FOR TEETH AND EARS



With DERMAVET LAVAPIROX® Vet Anti-Plaque C and DERMAVET LAVAPIROX® clean, you have two strong partners at your side to help you fight bacteria in the oral cavity and impurities in the inner ear.

Thanks to chlorhexidine, the DERMAVET LAVAPIROX® toothpaste fights plaque bacteria even without long exposure and brushing times, preventing tartar, inflammation and bad breath.

DERMAVET LAVAPIROX® Clean reliably frees the inner ear from impurities. Due to the high viscosity and the high proportion of squalane (23%), the ear is gently cleaned and the skin gently cared for. In addition, squalane dissolves stuck earwax and thus removes the nutritional basis for yeast fungi in particular.

! In case of bacterial infestation in the ear, we recommend the DERMAVET LAVAPIROX® Fluid for rinsing.

## COOL AND SOFT: FOR HOT DAYS AND DRY SKIN



### NEW AT DERMAVET LAVAPIROX®: THE COOLING GEL

Help your animal patients face the hot season with serenity. DERMAVET LAVAPIROX® Cool Gel helps regulate dogs' body temperature when summer shows its hottest side.

The formula with menthol and mint cools in a pleasant way, helping dogs to return more quickly to their usual shape in spite of the heat and intense physical strain, and calming the increased respiratory rate.

### NEW AT DERMAVET LAVAPIROX®: THE CREAM FOAM

Dry, chapped and cracked skin needs intensive care. DERMAVET LAVAPIROX® Keralyt forte helps stressed skin areas such as rough dog paws to become supple and elastic again.

Thanks to the formula with 15% urea, dry, brittle areas are optimally moisturized. This not only relieves unpleasant feelings of tension and itching, but can also prevent chapping and the associated risk of infection.





WE GIVE OUR BEST SO THAT YOU CAN GIVE YOUR BEST TOO

Even in times of increased hygiene and safety measures and heightened demand for protective equipment, we want to make sure that you have everything you need to work professionally and safely. Should an item be out of stock, we can offer you suitable replacements from other manufacturers - with the usual quality and speed of delivery.



AMPri Style by Med Comfort

Powder-free nitrile gloves, non-sterile

These gloves have a good grip, are skin-friendly and finger-textured.  
Glove length: 240 mm, layer thickness 0.08 mm.

Material: nitrile rubber

size	colour	pcs./pack. unit	Item No.
XS	apple	100 pieces	204584
S	apple	100 pieces	204170
M	apple	100 pieces	204171
L	apple	100 pieces	204172
XL	apple	100 pieces	204588

size	colour	pcs./pack. unit	Item No.
XS	hot chili	100 pieces	204549
S	hot chili	100 pieces	204550
M	hot chili	100 pieces	204551
L	hot chili	100 pieces	204552
XL	hot chili	100 pieces	204553

size	colour	pcs./pack. unit	Item No.
XS	prune	100 pieces	204554
S	prune	100 pieces	204555
M	prune	100 pieces	204556
L	prune	100 pieces	204557
XL	prune	100 pieces	204558

size	colour	pcs./pack. unit	Item No.
XS	grenadine	100 pieces	204559
S	grenadine	100 pieces	204560
M	grenadine	100 pieces	204561
L	grenadine	100 pieces	204562

size	colour	pcs./pack. unit	Item No.
XS	blueberry	100 pieces	204564
S	blueberry	100 pieces	204565
M	blueberry	100 pieces	204566
L	blueberry	100 pieces	204567
XL	blueberry	100 pieces	204568

size	colour	pcs./pack. unit	Item No.
XS	berry	100 pieces	204569
S	berry	100 pieces	204570
M	berry	100 pieces	204571
L	berry	100 pieces	204572
XL	berry	100 pieces	204573

size	colour	pcs./pack. unit	Item No.
XS	strawberry	100 pieces	204574
S	strawberry	100 pieces	204575
M	strawberry	100 pieces	204576
L	strawberry	100 pieces	204577

size	colour	pcs./pack. unit	Item No.
XS	grape	100 pieces	204579
S	grape	100 pieces	204580
M	grape	100 pieces	204581
L	grape	100 pieces	204582
XL	grape	100 pieces	204583

size	colour	pcs./pack. unit	Item No.
XS	orange	100 pieces	204589
S	orange	100 pieces	204590
M	orange	100 pieces	204591
L	orange	100 pieces	204592
XL	orange	100 pieces	204593

size	colour	pcs./pack. unit	Item No.
XS	lemon	100 pieces	204594
S	lemon	100 pieces	204595
M	lemon	100 pieces	204596
L	lemon	100 pieces	204597
XL	lemon	100 pieces	204598





PE veterinary gloves AMPri Med-Comfort

With extra long cuffs

The classic from AMPri for veterinarians and veterinary assistants as well as anyone who needs very long gauntlets. Made of polyethylene.  
Glove length: 900 mm. Packed in a bag.

Material: polyethylene

size	colour	pcs./pack. unit	Item No.
one size	orange	50 pieces	205085

Carton units: 40 packs

ADVANTAGES

- extra long cuff



AMPri Eco-Plus Latex

Powder-free, non-sterile latex gloves

AMPri Eco-Plus Latex is the all-rounder for all areas of application. The smooth surface gives this disposable glove a particularly good finger feel. It complies with the series of standards EN 420, EN 455 and EN 374. In addition, it is suitable for use with food in accordance with regulation (EC) 1935/2004 and complies with the specifications for personal protective equipment PPE category III in accordance with regulation (EU) 2016/425. Glove length: 240 mm

Material: latex

size	colour	pcs./pack. unit	Item No.
S	white	100 pieces	204916
M	white	100 pieces	204917
L	white	100 pieces	204918
XL	white	100 pieces	204919



AMPri Blue Eco-Plus Vinyl

Powder-free, non-sterile vinyl gloves

This disposable glove is powder-free, ambidextrous, non-sterile and produced according to EN 420, EN 455 and EN 374. Being free of latex proteins and accelerators, it is a good alternative to latex gloves, especially for people with type 1 latex allergies. Thanks to its special texturing, this glove is also characterized by a good grip. AMPri Blue Eco-Plus Vinyl is food safe according to regulation (EC) 1935/2004.

Material: vinyl

size	colour	pcs./pack. unit	Item No.
S	blue	100 pieces	204523
M	blue	100 pieces	204525
L	blue	100 pieces	204526
XL	blue	100 pieces	204527



Medi-Inn Nitril

Powder-free, non-sterile vinyl gloves

The microtextured fingertips of this glove provide a particularly good grip when handling dry and moist objects. The glove is also well suited for use in laboratories as PPE category 3 with protection against chemical hazards. As the glove is powder-free on the inside, there is no risk of contamination by powder during use. Glove length: 240 mm

Material: nitrile rubber

size	colour	pcs./pack. unit	Item No.
S	blue	200 pieces	205349
M	blue	200 pieces	205350
L	blue	200 pieces	205351
XL	blue	200 pieces	205352







Gentle Skin® Premium OP™

Powder-free surgical glove made of latex,  
sterile

This non-powdered surgical glove, sterile-packaged in pairs, has an anatomical fit that mimics the natural hand posture (left and right-handed), is micro-roughened on the fingers for sensitive tactile sensation and has a rolled edge.

Material: latex

size	colour	pcs./pack. unit	Item No.
Size 6	nature	50 pairs	1885
Size 6.5	nature	50 pairs	1886
Size 7	nature	50 pairs	1887
Size 7.5	nature	50 pairs	1888
Size 8	nature	50 pairs	1889
Size 8.5	nature	50 pairs	1890
Size 9	nature	50 pairs	1893

Carton units: 200 pairs

APPLICATION AREAS

- for all surgical areas with sterile requirements
- also suitable for double gloving



Reference OP®

Powdered surgical glove made of latex, sterile

This surgical glove, sterile-packaged in pairs, is lightly powdered, has an anatomical fit that mimics the natural hand posture (left and right-handed), is micro-roughened on the fingers for sensitive tactile sensation and has a rolled edge.

Material: latex

size	colour	pcs./pack. unit	Item No.
Size 6	nature	50 pairs	18675
Size 6.5	nature	50 pairs	18670
Size 7th	nature	50 pairs	18671
Size 7.5	nature	50 pairs	18672
Size 8	nature	50 pairs	18673
Size 8.5	nature	50 pairs	18674

Carton units: 300 pairs

APPLICATION AREAS

- for all surgical areas with sterile requirements
- also suitable for double-gloving



Suavel® Protec Plus

Highly filtering, disposable surgical mask  
with elastic ear loops, non-sterile

This 3-layer surgical mask is very skin-friendly thanks to the real surgical fleece and has a BFE of  $\geq 98\%$ .

size	colour	pcs./pack. unit	Item No.
one size	green	50 pieces	18728

Carton units: 900 pieces

APPLICATION AREAS

- for standard and long-term applications
- for all applications in the operating and practice area
- for the protection of users, visitors and patients



Suavel® Comfort Plus

Highly filtering, disposable surgical mask,  
to be tied, non-sterile


This 3-layer surgical mask is very skin-friendly thanks to the real surgical fleece, has extra long straps and a BFE of  $\geq 98\%$ .

size	colour	pcs./pack. unit	Item No.
one size	blue	50 pieces	18722

Carton units: 600 pieces

APPLICATION AREAS

- for standard and long-term applications
- for all applications in the operating and practice area
- for the protection of users, visitors and patients



Quality made in Germany:  
For a short video insight into  
Meditrade's mask production,  
scan the QR code on the right  
with your smartphone app.







### Suavel® Antifog

Highly filtering, disposable surgical mask for spectacle wearers, non-sterile

This 3-layer surgical mask is equipped with an effective spectacle-fog insulation in the form of a foam rubber insert. Thanks to the real surgical fleece, it is also very skin-friendly, has extra long straps and a BFE of  $\geq 98\%$ .

size	colour	pcs./pack. unit	Item No.
One size	blue	50 pieces	202826

Carton units: 600 pieces

#### APPLICATION AREAS

- for standard and long-term applications
- for all applications in the operating and practice area
- for the protection of users, visitors and patients



### Reusable gown "Krefeld"

Washable, 100 % cotton, 175 g / m<sup>2</sup>

This protective gown is made of 100 % high quality cotton (ISO certified, Ökotex standard 100 and without AZO dyes) and is washable at 95 degrees. The Krefeld model has a stand-up collar, knitted cuffs, a tie strap and loops at waist height for optimal adjustment. Very hard-wearing and long-lasting.

Material: 100 % cotton

size	colour	pcs./pack. unit	Item No.
XL - 130 cm x 140 cm	green	1 piece	205016
XL - 130 cm x 140 cm	blue	1 piece	205017

#### APPLICATION AREAS

- to protect the clothing of staff, visitors and patients from contamination
- unsterile

#### ADVANTAGES

- washable at 95 degrees



### BeeSana® PP protective sheath

Protective sheath made of PP fleece in the thickness 35 g / m<sup>2</sup>, non-sterile

The material is easy to wear and breathable as well as being liquid-repellent and durable – suitable for single use. Closure straps at on the neck and waist as well as tricot cuffs as sleeve endings guarantee good wearing comfort.

Material: polypropylene

size	colour	pcs./pack. unit	Item No.
140 cm x 140 cm	green	10 pieces	202836

Carton units: 100 pieces

#### APPLICATION AREAS

- to protect the clothing of staff, visitors and patients from contamination
- single use





BeeSana® PP gown

Protective gown made of PP fleece in the thickness 23 g / m<sup>2</sup>, non-sterile

Light, breathable material, tear-resistant and at the same time moisture-repellent – suitable for single use. Closure straps at the neck and waist as well as elastic rubber bands at the sleeves.

Material: polypropylene

size	colour	pcs./pack. unit	Item No.
110 cm x 140 cm	green	10 pieces	202829
120 cm x 140 cm	yellow	10 pieces	202831

Carton units: 100 pieces

APPLICATION AREAS

- to protect the clothing of staff, visitors and patients from contamination
- single use



BeeSana® SMS gown

Protective gown made of 3-layer, soft SMS composite material in thickness 40 g / m<sup>2</sup>, non-sterile

The combined layers guarantee a high barrier effect against liquids, germs and light splashes of chemicals – for single use. The gown is also opaque, lint-free, tear-resistant and breathable at the same time. Closure straps at the neck and waist as well as tricot cuffs as sleeve endings ensure a high level of comfort.

Material: polypropylene

size	colour	pcs./pack. unit	Item No.
115 cm x 150 cm	blue	10 pieces	202837

Carton units: 50 pieces

APPLICATION AREAS

- to protect the clothing of staff, visitors and patients from contamination
- single use



Suavel® Astrid

Universal beret hood

High-quality beret hood made of soft and light, breathable polypropylene non-woven fabric with a double-stitched elastic band. Good wearing comfort and pleasant fit.

Material: polypropylene

size	colour	pcs./pack. unit	Item No.
53 cm	blue	100 pieces	201789

APPLICATION AREAS

- for covering the hair



Surgical overshoes

Eco-Plus PE overshoes

Surgical overshoes made of polyethylene, elastic band made of bordered latex rubber, anti-slip effect.

Material: polyethylene

size	colour	pcs./pack. unit	Item No.
15 cm x 42 cm	blue	100 pieces	1870

Carton units: 5,000 pieces

APPLICATION AREAS

- for covering shoes



dermatest-  
garantie.de**VeraPanthen™****Hand & skin lotion with aloe vera & panthenol**

A rich care experience for every skin type with a delicate fresh scent. Thanks to the combination of innovative, highly concentrated and natural ingredients, VeraPanthen™ has a cell-renewing, anti-inflammatory and soothing effect. Its caring power is based on the optimal interaction of aloe vera gel in an exceptionally high concentration (25%), panthenol (provitamin B5) in a therapeutically effective concentration of 5% for inflammatory skin conditions and smoothing jojoba oil with its provitamin A and vitamin E.

Ingredients: aqua, aloe barbadensis leaf extract, iso-propyl-myristate, cetyl alcohol, panthenol, stearic acid, simmondsia chinensis (jojoba) seed oil, phenoxyethanol, piroctone olamine, benzoic acid, parfum, carbomer, dimethicone, polyacrylamide, C13-14 isoparaffin, laureth 7, triethanolamine, cetearyl alcohol, cetearyl glucoside, acid citric, BHT

size	colour	pcs./pack. unit	Item No.
180 ml	white	1 piece	202093

Carton units: 10 pieces, 84 pieces

**APPLICATION AREAS**

- is suitable for the care of closed wounds and inflammatory skin conditions thanks to the therapeutic concentration of panthenol
- for particularly dry skin and skin in need of care

**ADVANTAGES**

- absorbs quickly
- not sticky

**ROGG Skin Protection Cream****Perfume-free, water-in-oil protective cream**

This cream, with its skin-caring ingredients vitamin E, glycerine and beeswax, protects the skin from moisture loss, e.g. caused by frequent hand washing. The ROGG Skin Protection Cream restores the skin's natural oil-water protective film. Its protective effect has been confirmed according to the test methods in accordance with TRGS 401 / ABD and the recipe has been dermatologically tested.

Ingredients: aqua, paraffinum liquidum, glycerin, squalane, glyceryl oleate, cera alba, magnesium stearate, tocopheryl acetate, tocopherol, wheat amino acids, magnesium sulphate, hydrogenated palm glycerides citrate, citric acid, potassium sorbate, sodium benzoate

size	colour	pcs./pack. unit	Item No.
100 ml	white	1 piece	9717

Carton units: 12 pieces

**APPLICATION AREAS**

- skin protection for very dry skin areas

**ADVANTAGES**

- perfume free; the sensitive animal noses are not irritated
- restores the skin's natural protective oil-water film

**Fripa towel dispenser****For hand towels in the size 25 cm x 23 cm**

Robust dispenser made of metal, white-painted, for wall mounting.

size	colour	pcs./pack. unit	Item No.
27 cm x 13 cm x 35 cm	white	1 piece	185036

**APPLICATION AREAS**

- ideal for disposable towels

**Fripa disposable towels****2-ply, 25 cm x 23 cm**

Matching folded towels for dispensers with zigzag folds.

size	colour	pcs./pack. unit	Item No.
25 cm x 23 cm	white	3,000 pieces	202525

**Ecotex Plus disposable towels****2-ply, 25 cm x 23 cm**

Matching folded towels for dispensers with zigzag folds.

size	colour	pcs./pack. unit	Item No.
25 cm x 23 cm	green	3,750 pieces	20220

**ROGG dispenser box****Round, with visual control**

Transparent, elegant and practical dispenser box for e.g. pulp swabs or pulp-nonwoven compresses.

size	colour	pcs./pack. unit	Item No.
26 cm x 21 cm	transparent	1 piece	9804

Material: Plastic

**APPLICATION AREAS**

- for easy storage and removal of z. B. pulp swabs or pulp-nonwoven compresses



**Supplementary product:**  
ROGG instrument disinfection





Certoplast Prima insulating tape

Robust insulating tape

Robust black insulating tape for hoof and claw treatment, adhesive on both sides.

Material: viscose, rubber

size	colour	pcs./pack. unit	Item No.
45 mm x 25 m	black	1 roll	203641

APPLICATION AREAS

- for bandaging, applying and fixing protective bandages in the hoof and claw area



tesa textile adhesive tape

Robust textile adhesive tape

Robust, textile adhesive tape from tesa in four different colors.

Material: cellulose fabric, rubber

size	colour	pcs./pack. unit	Item No.
25 mm x 50 m	black	1 roll	203649
25 mm x 50 m	red	1 roll	203650
25 mm x 50 m	white	1 roll	203652

size	colour	pcs./pack. unit	Item No.
50 mm x 50 m	black	1 roll	203653
50 mm x 50 m	red	1 roll	203654
50 mm x 50 m	gray	1 roll	203656
50 mm x 50 m	white	1 roll	203655

APPLICATION AREAS

- for bandaging, applying and fixing protective bandages in the hoof and claw area



tesa masking tape

Water-repellent masking tape

Adhesive, solvent-free and writable fixing and masking tape from tesa. Made from 55 % renewable raw materials.

Material: rubber

size	colour	pcs./pack. unit	Item No.
19 mm x 50 m	light brown	1 roll	203657

APPLICATION AREAS

- for fixing bandages and dressings



tesa extra Power® Universal

Power film tape, duct tape

Extremely strong adhesive and fabric-reinforced film tape from tesa with a particularly fine structure for very high tear resistance. No fringing - even when simply tearing off by hand. Waterproof.

size	colour	pcs./pack. unit	Item No.
48 mm x 50 m	silver	1 roll	203659

APPLICATION AREAS

- for bandaging, applying and fixing protective bandages in the hoof and claw area





benartis silver tape professional

Particularly strong, robust foil fabric tape

Particularly strong, waterproof film tape, to be torn off by hand, does not fray.

Material:

size	colour	pcs./pack. unit	Item No.
5 cm x 50 m	silver	1 roll	7122

Carton units: 24 rolls

APPLICATION AREAS

- for bandaging, applying and fixing protective bandages in the hoof and claw area

benartis duct tape

Robust foil fabric tape

Robust textile adhesive tape in two different colors.

Material:

size	colour	pcs./pack. unit	Item No.
5 cm x 50 m	green	1 roll	7124
5 cm x 50 m	red	1 roll	7126

Carton units: 24 rolls

APPLICATION AREAS

- for bandaging, applying and fixing protective bandages in the hoof and claw area



mediware sterile surgical drapes

Without hole, 2-ply, waterproof

Comfortable and particularly user-friendly covers, for patient and inventory, made of particularly strong material. Completely liquid-impermeable, for the exact delimitation of the surgical field. Individually sterile in a peel pack and then packed in a dispenser box.

Material:

size	colour	pcs./pack. unit	Item No.
37 cm x 45 cm	blue	75 pieces	16799
90 cm x 150 cm	blue	25 pieces	200869
70 cm x 90 cm	blue	30 pieces	16792
45 cm x 75 cm	blue	50 pieces	16980

APPLICATION AREAS

- to cover the patient and inventory during all kinds of surgical interventions

ZARYS incision drape

Sterile, individually packed

Transparent, liquid-tight incision drape, with high breathability, which guarantees vapor permeability between the drape and the skin and prevents excessive moisture accumulation. It has good adhesive properties, is very thin and elastic and, therefore, adapts optimally to the shape of the body, thus guaranteeing a high degree of freedom of movement. Individually wrapped in paper foil.

Material: polyurethane, hypoallergenic acrylic adhesive

size	colour	pcs./pack. unit	Item No.
30 cm x 28 cm	transparent	20 piece	203675

APPLICATION AREAS

- to minimize the risk of infection and for the sterile sealing of the surgical field
- to protect sensitive areas during surgical interventions / operations





## Pulp / cellulose

Loose, unbleached, in layers

Crepe blanks in loose layers, made from recycled cellulose from environmentally friendly production.

Material: 100% recycled pulp

size	colour	pcs./pack. unit	Item No.
40 cm x 60 cm	unbleached	3 x 5 kg	1022

### APPLICATION AREAS

- as an underlay for non-sterile activities, as absorbent material or for cleaning and wiping



## Lohmann pulp

Loose, bleached, in layers

Cellulose blanks in loose layers, made from highly bleached cellulose.

Material: pulp

size	colour	pcs./pack. unit	Item No.
20 cm x 20 cm	bleached	3 x 5 kg	203715

### APPLICATION AREAS

- as an underlay for non-sterile activities, as absorbent material or for cleaning and wiping




## curea SafeCare

Strong hygiene underlay  
with SuperCore® absorption core

Strong hygiene pad with a skin-friendly surface and SuperCore® core for an uniform absorption and distribution of fluids. The patented technology from curea guarantees germ and odor-binding containment of liquids, high strength and stability - even under intensive use.

size	colour	pcs./pack. unit	Item No.
60 cm x 60 cm	white	112 pieces	204141
60 cm x 90 cm	white	70 pieces	204142



### APPLICATION AREAS

- as an underlay for non-sterile activities, as absorbent material as well as for post-operative care

### ADVANTAGES

- hypoallergenic and, therefore, suitable for sensitive animal noses thanks to the absence of adhesives and binding agents



Are you looking for  
sutures or cannulas?  
Please contact us!  
We would be pleased  
to order the items you  
need for you.





DERMAVET.  
LAVAPIROX®

Systematically  
against fungi,  
bacteria  
and itching

cleans,  
nourishes and  
eliminates  
odors



### DERMAVET LAVAPIROX® shampoo

For coat washing in case of bacterial  
and fungal infestation

Supportive treatment for bacterial, yeast infections, allergies, specific and non-specific itching. The antibacterial, fungicidal formula with salicylic acid fights germs, pathogens and fungi, typical dog odor and dandruff – all without drying out the skin, irritating mucous membranes or impairing wound healing.

Ingredients: Dehyquart Guar N, ethanol, sodium hydroxide, piroctone olamine, phenoxyethanol, Plantacare 818 UP (10%), polihexanide, salicylic acid, fragrance

size	pcs./pack. unit	Item No.
250 ml	1 bottle	205080

Carton unit: 9 pieces

#### APPLICATION AREAS

- for washing the coat of pets and farm animals such as cats, dogs and horses.
- weekly application for infestation with bacteria, yeast fungi, allergies, non-specific itching and dandruff
- for a lasting antibacterial, fungicidal treatment, we recommend between washes the application with DERMAVET LAVAPIROX® foam and DERMAVET LAVAPIROX® fluid

#### ADVANTAGES

- very well tolerated by skin and mucous membranes (also applicable in the area of ears, eyes and foreskin)
- no negative effect on wound healing
- without addition of artificial foaming agents, which could dry out the skin

1

coat wash once or twice a week for supportive therapy against fungi, bacteria and itching



### DERMAVET LAVAPIROX® foam

Coat cleaning without rinsing

DERMAVET LAVAPIROX® foam supports the treatment of bacterial and fungal infestations, as well as non-specific itching and allergies. It gently cleanses the coat and skin without additional moistening and rinsing. The antimicrobial active ingredients fight germs and incidentally also the typical dog smell.

Ingredients: Neutral alcohol, polyhexanide solution, Plantacare 818 UP, Octopirox, Phenoxtol

size	pcs./pack. unit	Item No.
200 ml	1 pumping bottle	205081

Carton unit: 9 pieces

#### APPLICATION AREAS

- Coat cleaning without rinsing for pets and farm animals such as cats, dogs and horses
- supportive therapy in case of infestation with bacteria, yeast fungi, allergies, unspecific itching and dandruff
- for water-free cleaning of stinky spots after rolling around in the dirt and to prevent the "dog smell" after swimming
- suitable for daily use
- as a complementary product to DERMAVET LAVAPIROX® shampoo for a lasting effect

#### ADVANTAGES

- use without water: especially suitable for water-shy animals and dry, irritated skin
- very well tolerated by skin and mucous membranes (can also be used in the area of ears, eyes and foreskin)
- no negative influence on wound healing
- fast drying of the treated areas

2

daily support for sustained bactericidal, fungicidal and keratolytic therapy



### DERMAVET LAVAPIROX® fluid

For spot treatment and rinsing

DERMAVET LAVAPIROX® fluid is the ideal complement to the coat wash with DERMAVET LAVAPIROX® shampoo. Since the antibacterial, fungicidal effect of the shampoo can only become active acutely during the weekly wash, we recommend an accompanying application with DERMAVET LAVAPIROX® fluid and DERMAVET LAVAPIROX® foam.

The fluid is particularly suitable for rinsing and spot treatment of infested areas, as well as for areas that are difficult to access, such as the foreskin and interdigital spaces.

Ingredients: 2-Phenoxyethanol, fragrance, ethanol, sodium hydroxide, piroctone olamine, Plantacare 818 UP, Polihexanide

size	pcs./pack. unit	Item No.
100 ml	1 bottle	205083

Carton unit: 9 pieces

#### APPLICATION AREAS

- complementary to the use of DERMAVET LAVAPIROX® Shampoo
- for spot treatment, rinsing and local disinfection of affected skin areas
- also for rinsing of foreskin and mucous membranes

#### ADVANTAGES

- funnel-shaped applicator facilitates precise application on affected skin areas
- very well tolerated by skin and mucous membranes (can also be used in the area of ears, eyes and foreskin)
- no negative influence on wound healing
- does not burn on irritated or sore skin
- not ototoxic

3

for local application, especially in hard to reach places and as a rinse





NEW!

DERMAVET LAVAPIROX® Keralyt forte

Intensive care for dry, chapped skin

Keralyt forte by DERMAVET LAVAPIROX® is a cream foam for intensive care of hard, brittle or chapped skin, such as on dog paws. The high content of urea (15%) makes dry, stressed skin soft and elastic again, which can prevent the development of itching, chapping and associated infections.

Ingredients: Aqua, Urea, Butane, Glycerin, C12-15 Alkyl Benzoate, C12-13 Alkyl Lactate, Cetearyl Alcohol, Persea Gratissima Oil, Propane, Cetearyl Glucoside, Hexyldecanol, Hexyldecyl Laurate, Sorbitol, Stearic Acid, Sodium Lauroyl Sarcosinate, Butyrospermum Parkii Butter, Isobutane, Hydroxyethyl Acrylate/Sodium Acryloyldimethyl Taurate Copolymer, Squalane, Lactic Acid, Sodium Lactate, Polysorbate 60, Sorbitan Isostearate

size	pcs./pack. unit	Item No.
125 ml	1 can with spray head	205310

Carton units: 24 pieces

APPLICATION AREAS

- care of hard, brittle, dry and cracked skin of pets and farm animals such as cats, dogs and horses

ADVANTAGES

- high percentage of urea
- without added fragrances or colouring agents



NEW!

DERMAVET LAVAPIROX® cool gel

Cooling gel for dogs with menthol and mint

DERMAVET LAVAPIROX® cool gel has been designed for the external use on dogs to help regulate their body temperature. In hot weather, the formula with menthol and mint can help to relieve the animal's discomfort, restlessness and increased breathing rate.

Ingredients: Alcohol, Aqua, Glycerin, Trometamol, Carbomer, Mentha Piperita Oil, Menthol

size	pcs./pack. unit	Item No.
250 ml	1 bottle	205082

Carton unit: 9 pieces

APPLICATION AREAS

- for external use on dogs
- for cooling at high outdoor temperatures

ADVANTAGES

- cooling effect and pleasant smell thanks to menthol and mint



DERMAVET LAVAPIROX® clean

For gentle cleaning of the ear canal

DERMAVET LAVAPIROX® clean can be used to gently clean the inner ear of pets and farm animals such as cats, dogs and horses. The high concentration of squalane (23%) not only gives DERMAVET LAVAPIROX® clean a viscous consistency that makes ear cleaning more pleasant for both patient and user, but also cares for the skin and makes it supple.

In case of bacterial infestation of the inner ear, the therapy can be supported by rinsing with the DERMAVET LAVAPIROX® fluid.

Ingredients: Paraffinum liquidum, 23 % Squalane, Isopropyl Myristate, Tocopherol

size	pcs./pack. unit	Item No.
100 ml	1 bottle	205083

Carton unit: 9 pieces

APPLICATION AREAS

- for cleaning dirty ear canals of pets and farm animals such as cats, dogs and horses
- eliminates unpleasant odors

ADVANTAGES

- high viscosity: the viscous formula makes DERMAVET LAVAPIROX® Clean easier to dose and distribute in the inner ear
- in addition, the high viscosity gives the product a more pleasant feeling in the animal patient, because it does not run into the ear immediately
- not ototoxic
- high concentration of squalane (23 %) nourishes the skin



DERMAVET LAVAPIROX® vet anti-plaque C

Antibacterial toothpaste for pets

Vet anti-plaque C by DERMAVET LAVAPIROX® is the only dental care for dogs and cats with the antibacterial active ingredient chlorhexidine. This helps fight or prevent bacteria that cause bad breath and plaque, as well as typical problems such as tartar, gingivitis and tooth loss. Thanks to the antibacterial formula DERMAVET LAVAPIROX® vet anti-plaque C is effective even with short brushing times.

Ingredients: Aqua, Hydrated Silica, Sorbitol, Calcium Carbonate, Hydroxyethylcellulose, Phenoxyethanol, Sodium Borohydride, Xylitol, Ethylhexylglycerin, Chlorhexidine Digluconate

size	pcs./pack. unit	Item No.
50 ml	1 tube	205084

Carton unit: 12 pieces

APPLICATION AREAS

- for daily teeth cleaning of small animals like dogs and cats
- antibacterial formula can prevent gingivitis, bad breath and tooth loss

ADVANTAGES

- effective against plaque bacteria thanks to chlorhexidine
- tasteless





Equital Arthro Star Plus Liquid

Building and maintaining of cartilage, joint casings, tendons and connective tissues

Equital Arthro Star Plus contains valuable ingredients to support the joint metabolism, to build up and maintain cartilage, joint casings, tendons and other connective tissues. It also contributes to the maintenance of the synovial fluid and thus supports the joint lubrication.

**Feeding recommendation:**  
Horse 30 ml / day  
Pony 15 ml / day  
Waiting time 96 h

⚠ **Feeding advice:** Devil's Claws should not be fed to pregnant mares.

📖 **Composition:** devil's claw brew (cut from devil's claw root) 40.00%, glucosamine hydrochloride 20.00%, methyl-sulfonylmethane (MSM) 5.00%, glycerine, willow bark brew (cut from willow bark), chondroitin sulfate, ground ginger root, hyaluronic acid

**Analytical components and contents:** crude protein 8.1%, crude fiber 0.0%, crude fats 0.0%, crude ash 0.0%, sodium 0.09%, moisture 74.7%

**Additives per liter:**  
Sensory additives: flavorings 2b: frankincense extract (boswellia serrata) 1000 mg

size	pcs./pack. unit	Item No.
1000 ml dosing bottle	1 piece	205022

ADVANTAGES

- with hyaluronic acid



Equital Broncho Horse Liquid

Assistance with nutritionally induced respiratory disorders

Broncho Horse Liquid contains valuable herbal extracts and black seed oil to support diet-related disorders of the respiratory tract.

**Feeding recommendation:**  
Horse 30 ml / day  
Pony 15 ml / day  
Waiting time 48 h

📖 **Composition:** 15% aqueous extract of herbs (fenugreek seeds whole, ribwort cut, echinacea herb cut, anise seeds rubbed, sage leaves cut, peppermint leaves cut, fennel seeds whole, yarrow herb cut, wormwood cut, nettle herb cut, thyme herb cut, watercress herb cut), 5% Licorice root juice, sucrose, 2% black seed oil, sodium caseinate

**Analytical components and contents:** crude protein 0.0%, crude fiber 0.0%, crude ash 0.0%, crude fats 1.4%, Humidity 93.6%, sodium 0.0%

size	pcs./pack. unit	Item No.
1000 ml dosing bottle	1 piece	205026



Equital Bone Star Powder

A balanced percentage of calcium and vitamin D3

By supplementing the feed with Equital Bone Star, deficiency symptoms can be compensated for in a meaningful way. With the usual hay and oat rations, the calcium content is much lower than the phosphorus content. This often causes a lack of calcium in horses.  
A balanced interaction of calcium and vitamin D3 can prevent:

- insufficient mineralization of the bone structure in young horses,
- a calcium deficiency in pregnant and lactating mares
- demineralization (osteomalacia and osteoporosis) in older horses.

Calcium is essential for the development, resilience and strength of bones and teeth.

Vitamin D3 promotes the absorption of calcium and phosphorus from the intestine and the reabsorption of calcium from the kidneys.

⚠ A longer feeding period with Bone Star is recommended for foals after weaning, as this lays the foundations for a healthy skeleton.

**Feeding recommendation for maintainance:**  
Small & large horse 50 g / day  
Foal / yearling 20-30 g / day  
Waiting time: none

**Feeding recommendation in case of deficiency symptoms:**  
Small & large horse 100 g / day  
Foal / yearling 50-60 g / day  
Waiting time: none

size	pcs./pack. unit	Item No.
3 kg bucket	1 piece	205025

ADVANTAGES

- can contribute to the mineralization of the bone
- particularly suitable in the growth phase
- relaxes the muscles and regenerates them in resting phases







Equital Hepa Forte Extra

Support of the liver function

The liver is an important metabolic organ and therefore plays a decisive role in the conversion of feed into performance-effective substances as well as in the detoxification and excretion of metabolic waste products. Equital Hepa Forte Extra can effectively support the liver in this process. The valuable ingredients from various herbs can intervene in diet-related metabolic problems. Organic sulfur improves the permeability of cell membranes. This improves the uptake of nutrients by the cells.

Feeding recommendation:

Horse 20 ml / day  
Pony 10 ml / day  
Waiting time: 48 h

**Composition:** fruit vinegar, 20% liquid herb extract (artichoke leaves cut, nettle leaves cut, ginseng root cut, milk thistle seeds rubbed), 5% methyl sulfonyl methane (MSM), elderberry juice, beetroot juice, fermented grain liquid, glucose


**Analytical components and contents:** crude protein 0.0%, crude fiber 0.0%, crude ash 0.0%, crude fats 0.0%, moisture 96.7 % sodium 0.06%

**Additives per liter:** sensory additives: flavorings 2b: 10,000 mg grapefruit seed extract (citrus aurantium)

size	pcs./pack. unit	Item No.
1000 ml dosing bottle	1 piece	205034

ADVANTAGES

- promotes the entire liver metabolism
- MSM increases the permeability of cell membranes



Equital Dermasan Horse Liquid

Improvement of the quality and strength of skin, hair and hoof horn

With Equital Dermasan Horse, skin and coat problems as well as hoof problems can be effectively avoided. A strong hoof is of great importance in horses. Problems in the hoof area often lead to lameness and severe pain. An optimal supply of biotin, zinc and the valuable omega-3 fatty acids EPA (0.8%) and DHA (0.5%) improves the quality and strength of skin, hair and hoof horn. Another important component for horn formation is sulfur. In addition to the sulfur-containing amino acid methionine, Equital Dermasan Horse contains MSM, a natural source of organic sulfur.

Feeding recommendation:

Horse 15 ml / day  
Pony 10 ml / day  
Waiting time: 48 h

**Composition:** vegetable glycerin, fish oil tasteless, methylsulfonylmethane (MSM) 5.0%, 1,2-propanediol, sodium caseinate

**Analytical components and contents:** crude protein 1.8%, crude fiber 0.0%, crude fats 3.7%, crude ash, 6.0%, methionine 0.18%, sodium 1.68%, moisture 73.9%

**Additives per liter:** nutritional additives: vitamin B1 (3a820) 175 mg, vitamin B2 as riboflavin 175 mg, vitamin B6 as pyridoxine hydrochloride (3a831) 250 mg, vitamin B12 as cyanocobalamin 750 mcg, niacinamide (3a315) 1,150 mg, biotin (3a880) 3,000,000 mcg, calcium D-pantothenate (3a841) 225 mg, folic acid (3a316) 50 mg, zinc as glycine zinc chelate, hydrate (3b607) 5,000 mg

size	pcs./pack. unit	Item No.
1000 ml dosing bottle	1 piece	205035



Equital Immun Aktiv Liquid

Supports the natural immune defense

Equital Immun Aktiv Liquid was specially designed to promote the natural immune defense. It contains vitamin E as well as vitamins of the B complex and essential minerals. The combination of these ingredients can reduce the concentration of toxins in the body. At the same time, Equital Immun Aktiv ensures that the electrolyte balance is not impaired, but rather remains stable. Overall, Equital Immun Aktiv promotes strong defenses for optimal well-being with its ideal combination of herbs, vitamins and minerals.

Feeding recommendation:

Horse 30 ml / day  
Pony 15 ml / day  
Waiting time: none

**Composition:** fruit vinegar, 1,2 propanediol, holdunderberry juice, beetroot juice, liquid fermented grain, liquid herbal extract from (artichoke leaves chopped, oregano herb rubbed, anise seeds rubbed, nettle herb chopped, thyme herb rubbed, milk thistle seeds rubbed, ribwort chopped), magnesium propionate, sodium diacetate, potassium chloride, magnesium chloride, calcium chloride

**Analytical components and contents:** crude fats 2.7%, crude ash 2.3%, crude protein 0.6%, crude fiber 0.0%, humidity 88.9%, sodium 0.63%

**Additives per liter:** vitamins, provitamins and chemically defined substances with a similar effect: (3a671) vitamin D3 32,000 IU, (3a700) vitamin E 10890 mg, niacinamide (3a315) 3,600 mg, vitamin B1 (3a820) 320 mg, vitamin B2 as riboflavin 480 mg, vitamin B6 as pyridoxole hydrochloride (3a831) 480 mg, vitamin B12 as cyanocobalamin 1800 mcg

**Compounds of trace elements:** Selenium as the hydroxy analog of selenomethionine (3b814) 4 mg; emulsifiers: lecithins (1c322) 15680 mg; antioxidants: butylated hydroxyanisole (BHA / E 320) 125 mg, propyl gallate (E 3101) 125 mg; flavorings 2b: 10,000 mg grapefruit seed extract (citrus aurantium)

size	pcs./pack. unit	Item No.
1000 ml dosing bottle	1 piece	205034

ADVANTAGES

- designed to strengthen the immune system
- promotes digestion
- regulates the metabolism in the case of sweet itch

## GENERAL TERMS AND CONDITIONS

### S1 Scope

(1) The following general terms and conditions (GTC) regulate the contractual relationship between ROGG Verbandstoffe GmbH & Co. KG and entrepreneurs, consumers, legal entities under public law or special funds under public law within the meaning of Section 310 (1) BGB. We only recognize conditions of the customer that contradict or deviate from our sales conditions if we expressly agree to their validity in writing.

(2) These conditions of sale also apply to all future transactions with the customer, insofar as they are legal transactions of a related nature.

(3) Conditions that deviate from our terms and conditions are only valid if we have expressly agreed to their validity. Deviations require a written agreement. If individual provisions of these conditions are changed or become ineffective, the remaining provisions remain in effect.

(4) The contract language is German.

### S2 Offer and conclusion of contract

(1) All our offers are non-binding with regard to price, quantity, delivery time and delivery options.

(2) Orders and online offers represent a binding offer to you to buy goods.

(3) If an order is to be regarded as an offer according to § 145 BGB, we can use this within one week by sending an order confirmation or by delivering the ordered goods to be taken.

(4) Orders are only legally effective when they have been confirmed by us in writing. We are not obliged to deliver due to the sending of price lists, circulars or orders received in response to general offers.

### S3 Documents provided

(1) For all documents provided to the customer in connection with the placing of the order – also in electronic form – such as price lists, calculations, drawings, etc., we reserve property rights and copyrights. These documents may not be made accessible to third parties unless we give the customer our express written consent. If we do not accept the offer of the customer within the period of section 2, these documents must be returned to us immediately.

(2) Telephone or electronic orders with details of the concluded contract (e.g., type of product, price, etc.) are saved by us. However, you cannot access your past orders over the Internet.

### S4 Prices and payment

(1) Unless otherwise agreed in writing, our sales prices are in € (plus VAT) and unless otherwise expressly agreed in writing, EXW (ex-works), excluding packaging.

(2) The payment of the purchase price is to be paid free of postage and expenses, in accordance with the payment conditions stated on the invoice, to one of our accounts listed on the invoices. We must have received the payment by the due dates. Discount deduction is only allowed upon written agreement and will be granted provided that all payment obligations resulting from previous shipments are met. Without a separate agreement, payment is to be made without any deduction.

(3) Unless otherwise agreed, the purchase price is to be paid within 30 days from the invoice date.. Interest on arrears will be charged at a rate of 5% for buyers and sellers and 8% for companies above the respective base rate pa (see Appendix 1). The assertion of a higher damage caused by default remains reserved. In the event that we assert higher damage caused by delay, the customer has the opportunity to prove to us that the asserted damage caused by delay did not occur at all or at least in a significantly lower amount.

(4) Interest is not charged on advance payments or payments on account. We reserve the right to cancel a loan, including such, within the payment deadlines according to the terms of payment at any time. The withholding or offsetting of payments on the part of the buyer due to any counterclaims is not permitted. Unless a fixed price agreement has been made, we reserve the right to make reasonable price changes due to changed wage, material and distribution costs for deliveries that are made 3 months or later after the conclusion of the contract.

### S5 Rights of retention

(1) The customer is only entitled to offset if his claims have been legally established or are undisputed. The purchaser is also entitled to offset against our claims if he asserts notices of defects or counterclaims from the same purchase contract. The customer is only authorized to exercise a right of retention if his counterclaim is based on the same contractual relationship.

### S6 Delivery time

(1) The start of the delivery time specified by us presupposes the timely and proper fulfilment of the purchaser's obligations. The exception of the unfulfilled contract remains reserved.

(2) Delivery times will be adhered to as best as possible but are non-binding. No claims, in particular no damage caused by delay, can be derived from late delivery.

(3) The consequences of force majeure of any kind, e.g., operational or traffic disruptions, fire damage, floods, deficiencies in work, energy, raw materials and auxiliary materials, strikes, lock-outs, disruptions in shipping, official measures and other circumstances for which we are not responsible, give us the right to stop deliveries without granting damages and without an obligation to make subsequent deliveries or to postpone delivery for the duration of the hindrance. Part deliveries are permitted. Call orders can only be carried out within the scope of the manufacturing possibilities.

(4) Deviations in size, weight, strength and quality are permissible within the scope of the customary market tolerances. Deviations from samples or earlier deliveries are avoided as far as technically feasible. The ordered quantities can be exceeded or undercut by up to 10%.

(5) If the customer is in default of acceptance or if he culpably breaches other obligations to cooperate, we are entitled to demand compensation for the damage we incur in this respect, including any additional expenses. We reserve the right to make further claims. If the above conditions are met, the risk of accidental loss or accidental deterioration of the purchased item is transferred to the customer at the point in time at which he is in default of acceptance or payment.

(6) Further legal claims and rights of the customer due to a delay in delivery remain unaffected.

### S7 Transfer of risk when shipping

(1) All shipments travel at the risk of the customer, regardless of whether carriage-paid delivery has been agreed. The seller reserves the right to choose the route and type of dispatch. Additional costs caused by special shipping requests by the buyer are at the buyer's expense, unless otherwise agreed. The same applies to increases in freight rates etc. occurring after the conclusion of the contract, unless carriage paid delivery has been agreed.

### S8 Retention of title

(1) We reserve title to the delivered item until all claims from the delivery contract have been paid in full. This also applies to all future deliveries, even if we do not always expressly refer to them. We are entitled to take back the purchased item if the customer behaves contrary to the contract.

(2) As long as ownership has not yet passed to him, the customer is obliged to treat the purchased item with care. As long as ownership has not yet passed, the customer must notify us immediately in writing if the delivered item is seized or exposed to other interventions by third parties. If the third party is unable to reimburse us for the judicial and extrajudicial costs of a lawsuit in accordance with Section 771 of the German Code of Civil Procedure, the customer is liable for the loss we incur.

(3) The customer is entitled to process or sell goods in the ordinary course of business even before the purchase price has been paid. The customer hereby assigns to us the claims against the customer from the resale of the goods subject to retention of title in the amount of the final invoice amount agreed with us (including VAT). This assignment applies regardless of whether the purchased item has been resold without or after processing.

(4) The buyer must immediately notify us of any seizure or any other impairment of our property rights, confirm the property rights to both third parties and us in writing and support us in safeguarding our rights in every way. If the buyer sells the goods delivered by us, regardless of the condition, he assigns the claims arising from the sale to his customers with all ancillary rights to us with the conclusion of the contract until the complete repayment of all our claims, namely in the amount of the invoice value of our reserved goods contained in the sold items. At

our request, the buyer is obliged to notify his customers of the assignment and to give us the information necessary to assert our rights and to hand over the necessary documents.

(5) The customer remains authorized to collect the claim even after the assignment. Our authorization to collect the claim ourselves remains unaffected. However, we will not collect the claim as long as the customer fulfils his payment obligations from the proceeds received, is not in default of payment and, in particular, no application has been made to open insolvency proceedings or payments have been suspended.

(6) If, after the conclusion of the contract, circumstances become known that give rise to concerns about the creditworthiness of the buyer or if the buyer defaults on the fulfilment of his obligations towards us, we are entitled to immediate payment or, for our security, the surrender of the goods and to request advance payment or security for future goods to be delivered. The filing of bankruptcy, the taking of the oath of disclosure, a change of the company owner in connection with payment difficulties release us from the fulfilment of any ongoing orders, whereby the decision on the execution of the orders lies with us.

### S9 Warranty and notification of defects as well as recourse / manufacturer recourse

(1) Guarantee to consumers

a) Insofar as the information contained in our brochures, advertisements and other offer documents has not been expressly designated as binding by us, the illustrations or drawings contained therein are only approximate.

b) Claims for defects become statute-barred 12 months after delivery of the goods we have delivered to our purchaser in the case of used goods. The statutory limitation period applies to claims for damages in the event of intent and gross negligence as well as injury to life, body and health based on an intentional or negligent breach of duty on the part of the user. Insofar as the law prescribes longer periods in accordance with Section 438 (1) No. 2 BGB (buildings and items for buildings), Section 445 b BGB (right of recourse) and Section 634a (1) BGB (construction defects), these periods shall apply. Prior to returning the goods our permit is to be requested.

c) If the delivered item does not have the quality agreed between the customer and us, or if it is not suitable for the use stipulated in our contract or for general use, or if it does not have the properties that the customer could expect according to our public statements, so we are obliged to provide supplementary performance. This does not apply if we are entitled to refuse supplementary performance due to legal regulations.

d) The customer initially has the choice of whether the supplementary performance should be carried out through repair or replacement. However, we are entitled to refuse the type of supplementary performance chosen by the customer if it is only possible at disproportionate costs and the other type of supplementary performance does not result in significant disadvantages for the customer. Reduction of the purchase price or withdrawal from the contract by the customer are excluded during subsequent performance. A subsequent improvement is deemed to have failed with the unsuccessful second attempt, unless something else results in particular from the nature of the item or the defect or the other circumstances. If the supplementary performance has failed or if we have refused the supplementary performance altogether, the customer can, at his option, request a reduction in the purchase price (reduction) or withdraw from the contract.

e) The customer can only assert claims for damages under the following conditions due to the defect if the supplementary performance has failed or we have refused the supplementary performance. The right of the customer to assert further claims for damages under the following conditions remains unaffected.

(2) Guarantee towards entrepreneurs

a) The purchaser's warranty rights presuppose that he has properly complied with his inspection and complaint obligations under Section 377 of the German Commercial Code (HGB).

b) Claims for defects become statute-barred 12 months after delivery of the goods delivered by us to our customer. The statutory limitation period applies to claims for damages in the event of intent and gross negligence as well as injury to life, body and health based on an intentional or negligent breach of duty on the part of the user. Insofar as the law prescribes longer periods in accordance with Section 438 (1) No. 2 BGB (buildings and items for buildings), Section 445 b BGB (right of recourse) and Section 634a (1) BGB (construction defects), these periods shall apply. Prior to returning the goods our permit is to be requested.

c) If, despite all due care, the delivered goods show a defect that already existed at the time of the transfer of risk, we will either repair the goods or deliver replacement goods, subject to timely notification of defects. We must always be given the opportunity to provide supplementary performance within a reasonable period of time. Recourse claims remain unaffected by the above regulation without restriction.

d) If the supplementary performance fails, the customer can – without prejudice to any claims for damages – withdraw from the contract or reduce the remuneration.

e) Claims for defects do not exist in the case of only insignificant deviations from the agreed quality, in the case of only insignificant impairment of usability, in the case of natural wear and tear as in the case of damage that after the transfer of risk as a result of incorrect or negligent treatment, excessive use, unsuitable operating resources, defective construction work, unsuitable building ground or due to special external influences that are not required by the contract. If improper repair work or changes are carried out by the customer or a third party, no claims for defects exist for these or the consequences arising from them.

f) Claims of the customer due to the expenses necessary for the purpose of the supplementary performance, in particular transport, travel, labour and material costs, are excluded if the expenses increase because the goods delivered by us are subsequently transferred to a different location than the purchaser's establishment unless the shipment corresponds to its intended use.

g) The purchaser's right of recourse against us only exists insofar as the purchaser has not made any agreements with his customer that go beyond the legally mandatory claims for defects. Concerning the extent of the customer's recourse entitlements against the supplier, paragraph 6 shall be applied accordingly.

### S10 Limitation of Liability

(1) We exclude liability for slightly negligent breaches of duty.

(2) Without prejudice to the above regulations and the following limitations of liability, we are fully liable for damage to life, body and health based on a negligent or wilful breach of duty by our legal representatives or our vicarious agents, as well as for damage covered by liability under the Product Liability Act , as well as for all damages based on wilful or grossly negligent breaches of contract as well as malice, our legal representatives or our vicarious agents. Insofar as we have given a quality and / or durability guarantee with regard to the goods or parts thereof, we are also liable under this guarantee. For damage that is based on the lack of the guaranteed quality or durability, but does not occur directly on the goods, we are only liable if the risk of such damage is clearly covered by the quality and durability guarantee.

(3) We are also liable for damage caused by simple negligence, insofar as this negligence relates to the breach of contractual obligations, compliance with which is of particular importance for achieving the purpose of the contract (cardinal obligations). However, we are only liable insofar as the damage is typically associated with the contract and is foreseeable. Otherwise, we are not liable for simple negligent breaches of secondary obligations that are not essential to the contract. The limitations of liability contained in sentences 1 – 3 also apply insofar as the liability for the legal representatives, executives and other vicarious agents is affected.

(4) Any further liability is excluded regardless of the legal nature of the asserted claim. Insofar as our liability is excluded or limited, this also applies to the personal liability of our employees, workers, employees, representatives and vicarious agents.

### S11 Miscellaneous

(1) This contract and all legal relationships between the parties are subject to the law of the Federal Republic of Germany to the exclusion of the UN Sales Convention (CISG).

(2) The place of performance and exclusive place of jurisdiction and for all disputes arising from this contract is our place of business, unless otherwise stated in the order confirmation.

(3) All agreements made between the parties for the purpose of executing this contract are set out in writing in this contract.

## RIGHT OF WITHDRAWAL

## RIGHT OF WITHDRAWAL

You have the right to withdraw from this contract within fourteen days without giving any reason.

The cancellation period is fourteen days from the day on which you or a third party named by you, who is not the carrier, has taken possession of the goods.

In order to exercise your right of withdrawal, you must inform us (ROGG Verbandstoffe GmbH & Co. KG, Oberndorfer Str. 7, 85777 Fahrenzhausen, [bestellung@roggvet.de](mailto:bestellung@roggvet.de), telephone: 49 8133 9320-0, fax: 49 8133 9320-25) by means of a clear Declaration (e.g. a letter sent by post, fax or email) of your decision to withdraw from this contract. You can use the attached model withdrawal form for this purpose, but this is not mandatory.

In order to meet the withdrawal deadline, it is sufficient for you to send your notification of exercising your right of withdrawal before the withdrawal period has expired.

## CONSEQUENCES OF WITHDRAWAL

If you withdraw from this contract, we will have given you all payments that we have received from you, including delivery costs (with the exception of the additional costs that result from choosing a different type of delivery than the cheapest standard delivery offered by us have), to be repaid immediately and at the latest within fourteen days from the day on which we received notification of your cancellation of this contract. For this repayment, we will use the same means of payment that you used for the original transaction, unless something else was expressly agreed with you; In no case will you be charged any fees for this repayment. We can refuse repayment until we have received the goods back or until you have provided evidence that you have sent the goods back, whichever is earlier. You must return or hand over the goods to us immediately and in any case no later than fourteen days from the date on which you informed us of the cancellation of this contract. The deadline is met if you send the goods before the period of fourteen days has expired.

You bear the direct costs of returning the goods.

You only have to pay for any loss in value of the goods if this loss in value is due to handling of the goods that is not necessary to check the nature, properties and functionality of the goods.



## WITHDRAWAL FORM

(If you want to cancel the contract, please fill out this form and send it back.)

To  
**ROGGvet**  
Oberndorfer Str. 7  
85777 Fahrenzhausen  
Germany

Phone: +49 (0)8133 9320-10  
Fax: +49 (0)8133 9320-25

Mail: [sales@roggvet.de](mailto:sales@roggvet.de)

I / we (\*) hereby revoke the purchase contract concluded by me / us (\*)  
of the following goods (\*) / the provision of the following service (\*)

\_\_\_\_\_  
\_\_\_\_\_

Ordered on \_\_\_\_\_ (\*)/received at \_\_\_\_\_ (\*)

Name of the consumer(s) \_\_\_\_\_

Address of the consumer(s) \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

\_\_\_\_\_  
Date Signature of the consumer (s) (only if this is communicated on paper)

(\*) Delete where inapplicable

## ROGGvet

Oberndorfer Str. 7  
85777 Fahrenzhausen  
Germany

## COPY TEMPLATE

**FAX: +49 (0)8133/9320-25**  
**[sales@roggvet.de](mailto:sales@roggvet.de)**

## ORDER FORM

We order according to the delivery and payment conditions known to us:

Article no.	Quantity	Article name	Packing unit / piece	Total €

Total:

Customer number.: .....

Date: .....

Signature: .....

Signature:

All prices plus VAT – free delivery from a net order value of 50 €.  
The general terms and conditions of ROGGvet, represented  
by ROGG Verbandstoffe & Co. KG, apply.



## Sustainable print products.



Our print products are:

- produced and shipped climate neutrally
- made with certified recycled paper
- printed with vegan inks made from renewable, raw plant materials



## Strong Brands.

